


# Fastorb

Surgical Suture, Polyglycolicacid (PGA)  
Fast Absorbable, Synthetic, Braided, Undyed

**Surgical Sutures**  
**Absorbable**



**BOZ**<sup>®</sup>

## Definition

Synthetic absorbable, coated, short term wound support and short term mass absorption.

## Material

Polyglycolic Acid (PGA)

## Coating

Calcium Stearate and polycaprolactone

## Structure

Braided

## Colour

Beige (undyed)

## Size

EP : 0.4 to 8

USP : 8/0 to 6

## Tensile Strength Retention

7 Days 75 %

## Mass Absorption

42 days

## Sterilization Method

Gamma Irradiation

## Main Indication

Gynaecology (episiotomies)

Urology

Ophthalmic surgery (conjunctive suturing)

Odontology (oral mucosa)























Pediatric surgery

Skin closer - intracutaneous, subcutaneous.

Ligatures

















# Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color	
<b>12 mm</b>						
<b>3/8</b> Reverse Cutting <b>330 Micron</b>	 	<b>FS1545</b>	<b>6/0</b> 0,7	75	12	Undyed
<b>13 mm</b>						
<b>1/2</b> Round Body, Taper Point <b>330 Micron</b>	 	<b>FS1563</b>	<b>6/0</b> 0,7	75	12	Undyed
<b>1/2</b> Round Body, Taper Point <b>380 Micron</b>	 	<b>FS1886</b>	<b>5/0</b> 1,0	75	12	Undyed
<b>3/8</b> Reverse Cutting <b>330 Micron</b>	 	<b>FS5912</b>	<b>6/0</b> 0,7	75	12	Undyed
<b>3/8</b> Reverse Cutting <b>380 Micron</b>	 	<b>FS1929</b>	<b>5/0</b> 1,0	75	12	Undyed
<b>3/8</b> Reverse Cutting <b>Premium 380 Micron</b>	 	<b>FS1941</b>	<b>5/0</b> 1,0	45	12	Undyed
<b>16 mm</b>						
<b>1/2</b> Round Body, Taper Point	 	<b>FS1989</b>	<b>5/0</b> 1,0	75	12	Undyed
		<b>FS2503</b>	<b>4/0</b> 1,5	75	12	Undyed
<b>3/8</b> Reverse Cutting	 	<b>FS1673</b>	<b>6/0</b> 0,7	75	12	Undyed
		<b>FS2019</b>	<b>5/0</b> 1,0	75	12	Undyed
		<b>FS2539</b>	<b>4/0</b> 1,5	75	12	Undyed
<b>3/8</b> Inside Cutting	 	<b>FS2557</b>	<b>4/0</b> 1,5	75	12	Undyed
		<b>FS3038</b>	<b>3/0</b> 2,0	75	12	Undyed
<b>19 mm</b>						
<b>3/8</b> Reverse Cutting	 	<b>FS2639</b>	<b>4/0</b> 1,5	75	12	Undyed
		<b>FS3108</b>	<b>3/0</b> 2,0	75	12	Undyed
<b>20 mm</b>						
<b>1/2</b> Round Body, Taper Point	 	<b>FS2228</b>	<b>5/0</b> 1,0	75	12	Undyed






# Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>20 mm</b>							
<b>1/2</b> Round Body, Taper Point			<b>FS2683</b>	<b>4/0</b> 1,5	75	12	Undyed
			<b>FS3141</b>	<b>3/0</b> 2,0	75	12	Undyed
<b>3/8</b> Reverse Cutting			<b>FS2250</b>	<b>5/0</b> 1,0	75	12	Undyed
<b>22 mm</b>							
<b>1/2</b> Round Body, Taper Point			<b>FS2745</b>	<b>4/0</b> 1,5	75	12	Undyed
<b>25 mm</b>							
<b>1/2</b> Round Body, Taper Point			<b>FS2787</b>	<b>4/0</b> 1,5	75	12	Undyed
			<b>FS3251</b>	<b>3/0</b> 2,0	75	12	Undyed
<b>3/8</b> Reverse Cutting			<b>FS3273</b>	<b>3/0</b> 2,0	75	12	Undyed
			<b>FS3756</b>	<b>2/0</b> 3,0	75	12	Undyed
<b>26 mm</b>							
<b>1/2</b> Round Body, Taper Point			<b>FS3785</b>	<b>2/0</b> 3,0	75	12	Undyed
			<b>FS4214</b>	<b>0</b> 3,5	75	12	Undyed
<b>3/8</b> Reverse Cutting			<b>FS3821</b>	<b>2/0</b> 3,0	75	12	Undyed
			<b>FS4236</b>	<b>0</b> 3,5	75	12	Undyed
			<b>FS3345</b>	<b>3/0</b> 2,0	75	12	Undyed
<b>30 mm</b>							
<b>1/2</b> Round Body, Taper Point			<b>FS3368</b>	<b>3/0</b> 2,0	75	12	Undyed
			<b>FS3857</b>	<b>2/0</b> 3,0	75	12	Undyed
			<b>FS4243</b>	<b>0</b> 3,5	75	12	Undyed
			<b>FS4607</b>	<b>1</b> 4,0	75	12	Undyed







# Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>30 mm</b>					
<b>3/8</b> Reverse Cutting 	<b>FS2885</b>	<b>4/0</b> 1,5	75	12	Undyed
	<b>FS3396</b>	<b>3/0</b> 2,0	75	12	Undyed
	<b>FS4628</b>	<b>1</b> 4,0	75	12	Undyed
<b>35 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS3915</b>	<b>2/0</b> 3,0	75	12	Undyed
	<b>FS4299</b>	<b>0</b> 3,5	75	12	Undyed
	<b>FS4305</b>	<b>0</b> 3,5	90	12	Undyed
	<b>FS4643</b>	<b>1</b> 4,0	75	12	Undyed
	<b>FS4648</b>	<b>1</b> 4,0	90	12	Undyed
<b>3/8</b> Reverse Cutting 	<b>FS3422</b>	<b>3/0</b> 2,0	75	12	Undyed
	<b>FS3932</b>	<b>2/0</b> 3,0	75	12	Undyed
<b>37 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS3954</b>	<b>2/0</b> 3,0	90	12	Undyed
<b>40 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS3981</b>	<b>2/0</b> 3,0	75	12	Undyed
	<b>FS3988</b>	<b>2/0</b> 3,0	90	12	Undyed
	<b>FS4367</b>	<b>0</b> 3,5	75	12	Undyed
	<b>FS4374</b>	<b>0</b> 3,5	90	12	Undyed

# Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>40 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS4712</b>	<b>1</b> 4,0	75	12	Undyed
	<b>FS4718</b>	<b>1</b> 4,0	90	12	Undyed
	<b>FS5001</b>	<b>2</b> 5,0	75	12	Undyed
<b>3/8</b> Reverse Cutting 	<b>FS4005</b>	<b>2/0</b> 3,0	75	12	Undyed
<b>45 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS4421</b>	<b>0</b> 3,5	75	12	Undyed
	<b>FS4771</b>	<b>1</b> 4,0	75	12	Undyed
	<b>FS5034</b>	<b>2</b> 5,0	75	12	Undyed
<b>Straight</b> Reverse Cutting 	<b>FS2911</b>	<b>4/0</b> 1,5	75	12	Undyed
<b>50 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS4037</b>	<b>2/0</b> 3,0	75	12	Undyed
	<b>FS4453</b>	<b>0</b> 3,5	75	12	Undyed
	<b>FS4811</b>	<b>1</b> 4,0	75	12	Undyed
	<b>FS5065</b>	<b>2</b> 5,0	75	12	Undyed
<b>65 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>FS4864</b>	<b>1</b> 4,0	100	12	Undyed

# Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>150 cm</b>					
<b>FASTSORB</b> Non-Needle Suture	<b>FS1373</b>	<b>6/0</b> 0,7	150	12	Undyed
	<b>FS2325</b>	<b>4/0</b> 1,5	150	12	Undyed
	<b>FS2934</b>	<b>3/0</b> 2,0	150	12	Undyed
	<b>FS3510</b>	<b>2/0</b> 3,0	150	12	Undyed
	<b>FS4113</b>	<b>0</b> 3,5	150	12	Undyed
	<b>FS4531</b>	<b>1</b> 4,0	150	12	Undyed
	<b>FS4907</b>	<b>2</b> 5,0	150	12	Undyed