

Tende Babyblue

LED Phototherapy Device for Neonates

Tende Babyblue is Turkey's first intensive LED phototherapy unit.


This device has powerful LED lights and it is used for the treatment of neonates with high bilirubin levels.













Top of the incubators



Angle adjustment (360°)

FEATURES

LIGHT SOURCE

Light Source Blue LED (24 pcs, Super Bright LED)

Wavelength 460nm (± 2%)

Light Intensity (Standard) Max. $60 \mu w/cm^2/nm (\pm 10\%)$ Min. $15 \mu w/cm^2/nm (\pm 10\%)$ (From 40 cm) Intensive Light Intensity (Optional) Max. $120 \mu w/cm^2/nm (\pm 10\%)$ Min. $30 \mu w/cm^2/nm (\pm 10\%)$ (From 40 cm)

tensive Light Intensity (Optional) wax. 120 μw/cm/nm (±10%) (From 40 cm

Effective Area 45 cm x 25 cm (From 40 cm)

LED Life Time > 70,000 Hours (Check the user manual for technical details)

ELECTRICAL SPECIFICATIONS

Operating Voltage ~110 / 240 V 60 /50 Hz (MDD Class/Safety)

Power Consumption 80 VA, Max 45 W

ENVIRONMENTAL REQUIREMENTS

Operating Temperature 15°C to 40°C
Operatin R. Humidity <90%

Storage Temperature ~ -20°C to +50°C

Storage R. Humidity < 90%

GENERAL

LCD Graphic Display > 5 inches

Dimensions L W H (Therapy Module) \sim 47 cm x 27.5 cm x 8 cm Wheels 2/3 with Brakes Height Adjustment 110 cm - 154 cm (\pm 2%)

Height Adjustment 110 cm - 154 cm (± Tilt Adjustment (with Stand Unit) Yes

Wheel Height ~ 6.3 cm

Weight $^{\sim}$ 4 Kg (11.2 Kg with Stand Unit) \pm 1%

APPLICATION

Treatment Time Adjustment Yes (Hours and Minutes hhh:mm)

Visual Alarms
Yes
Audible Alarms
Yes
Light Intensity
Date & Time Display
Yes

LED Life Counter Yes (Individual Level and Total LED Counters)

Baby Skin Temperature Yes (Optional)

MEDICAL DEVICES DIRECTIVES (MDD)

Manufacturer MDD Conformity

Class IIb - Production Quality Assurance
Equipment Conformity

EN60601-1, EN60601-1-2, EN60601-2-50



TENDE ELEKTRONİK YAZILIM MÜH. İLETİŞİM MAK. SAN. VE TİC. LTD. ŞTİ.

Ahi Evran OSB Mahallesi, Akhun Caddesi No: 5/1, 06935 Sincan ANKARA TÜRKİYE T: +90 312 397 0733 F: +90 312 397 1787 www.tende.com.tr info@tende.com.tr