

BOJIN MEDICAL www.bojin-medical.com



System 4200 Multi Functional Surgical Power Tools

For Joint surgery | Trauma surgery | Neurosurgery | Thoracic surgery

ATTACHMENTS	Model	Indication	Description	Speed	Cannulation
	BJ4201	Trauma- prosthesis	Sagittal saw use Bojin type blades	0 - 15000 cycles/min	-
	BJ4202	Trauma	0.5 - 7.4mm 1/4"	0 - 1100 rpm	-
	BJ4202sao	Trauma	A.O. quick coupling	0 - 1100 rpm	
	ВЈ4203В	Trauma	Zimmer reamer for nail surgery	0 - 1100 rpm 3.0 Nm	4.2mm
	BJ4203D	Trauma	A.O. reamer for nail surgery	0 - 1100 rpm 3.0 Nm	4.2mm
	BJ4204	Neuro	Self-stop cranial drill	0 - 1100 rpm	1
	BJ4205	Neuro	Cranial bur use together with cranial drill	0 - 40000 rpm	-
	BJ4207B	Prosthesis	Zimmer high torque reaming drill	0 - 300 rpm 8.0 Nm	4.2mm
	BJ4207D	Prosthesis	A.O. high torque reaming drill	0 - 300 rpm 8.0 Nm	4.2mm
=(BJ4209	Prosthesis	Reciprocating saw	0 - 15000 cycles/min	-
	BJ4213	Trauma	Wire driver A Wire driver B	0 - 1100 rpm 0 - 1100 rpm	0.8 - 2.0mm 2.0 - 3.0mm
	ACCESSOR	RIES			
					Y
	BJ4001 Battery	BJ4002 Charger	BJ4003 Cover	BJ4004 Aseptic transfer kit	BJ1005 Key

Auction no. 21075923 / ocds-b3wdp1-MD-1678893907083 Technical specifications

Electric Sagital Saw (Lot 15 & 16)

Auction specifications	BOJIN Specifications			
Hand Piece - BJ4200				
Description of Function: The Battery Operated Saw Attachment for sawing use in general traur	d Drill system is a battery driven power system with a Sagital matology.			
Electric Hand piece: Yes	Electric Hand Piece (Battery operated)			
Power: ≥ 100 W	Power: ≥ 100 W			
Weight: ≤ 950 g	Weight: ≤ 950 g			
Material: Stainless Steel	Motor body: Stainless Steel			
Sterilization: Autoclavable 135°	Sterilization method: Autoclavable 135°			
Speed: ≥ 900 rpm	Adjustable speed: 0 - 1100 rpm			
Noice: ≤ 75 db	Noice: ≤ 75 db			
Operating time: ≥ 15 min	Operating time: ≥ 30 min			
Blocking mechanism: Yes	Battery & Attachment blocking mechanism			
Attachment coupling: easily removable	Quick Coupling Attachments			
Increasing temperature while working: ≤ 20°C	Increasing temperature while working: ≤ 20°C			
Sterilizable cover: Yes	Hand Piece cover, Battery cover & Aseptic battery transfer kit: Autoclavable 135°			
	Forward and Reverse switch with safe mode			
	Pistol grip handpiece, ergonormic design			
Battery - BJ4001	是一种的一种。 第二十二十二十二十二十二十二十二十二十二十二十二十二十二十二十二十二十二十二十			
Quantitus 2 man, appilly removable	2 pcs.			
Quantity: 2 pcs., easily removable	Insert/detachment method: telescopic, easily			
Run time: ≥ 100 min	Run time with fully charged battery: ≥ 30 min			
	Ni-Mh battery, 14.4V/1800mAh			
	Non-memory			
	Non-autoclavaible			
	Shield to protect battery from the housing			
	Blocking mechanism: Yes			
Battery Charger - BJ1002-4				
Simultaneosly charged batteries: min. 2 batteries	4 batteries			
Battery charge status: Yes	charge status: Yes Light indicator			
Main power: 220V / 50Hz	Input Voltage: 110-250V, 50/60Hz.			
·	The output voltage: 17.5*4 (this is for four charger).			
	Output Current: 700mAh			

Sagittal Saw Attachment - BJ4201			
Sawing Speed: 0 - 20000 cpm	Adjustable speed: 0 - 15000 cycles/min		
Saw blade coupling: without key	Tool less blade mounting		
	360° rotating: Up, down, left and right		
	Sterilization method: Autoclavable 135°		
	Material: Stainless Steel		
	Blade arc of excursion: 5 degree		
Accessories			
Power cable: Yes	European power cable (battery charger)		
Container with fixing accessories during sterilization	Code: BJXDH4200. Container for sterilization, transportation and storage,		
Case for storage and transportation	with fixing accessories during sterilization		
Saw blade wide 10mm	Code: BJ90101.0A. Width: 10mm; Length: 90mm; Thickness: 1.0mm. Quantity: 5 pcs.		
Saw blade wide 15-20mm	Code: BJ90181.0A. Width: 18mm; Length: 90mm; Thickness: 1.0mm. Quantity: 10 pcs.		
Saw blade wide 25-30mm	Code: BJ90251.0A. Width: 25mm; Length: 90mm; Thickness: 1.0mm. Quantity: 10 pcs.		
Other requirements			
Warranty: 2 years	Handpiece, Attachments & Battery Charger: 2 years warranty		

Date: 12.04.2023

Name:

sition: Departu

3

BJ1002-4 Battery Charger (four chargers)



Fit for BJ1001, BJ4001, BJ6001 battery.

Input Voltage: 110-250V, 50/60Hz.

The output voltage: 17.5*4 (this is for four charger).

Output Current: 700mAh.

Battery charge status: Light indicator.

Size: 26(L)*26(W)*11cm(H)

Net weight: 1.4 kgs

Package size: 40(L)*28(W)*18(H)

GW: 2kgs

The manufacturers listed in this document are not affiliated with Bojin and have not endorsed the Bojin products listed in this catalogue.

This documentation is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part.

The information and illustrations included in this catalogue are correct and true to the best of our knowledge. However, we cannot be regarded as binding warranty for quality. The user is responsible for checking the products prior to use whether they are suited for the intended purpose. We can only be held liable for these instruments pursuant to our terms of delivery and payment. Subject to






Shanghai Bojin Medical Instrument Co.,Ltd

Company address: Room 606, No.3388, Congnesin Road, Shanghai, China Factory address: I-D, 2-A, No.128, Longpan Road, Jioding District, Shanghai, China Tel: 186-21-64538678.64538077.6652801

Fax:+86-21-66527863

Website: www.bajin-medical.com Email: manager@bojin-medical.com



ORTHOPEDIC CONSUMABLES

OVERVIEW









