

TECHNICAL SHEET

BONE INJECTION NEEDLES BEVELED TIPS — PLASTIC

Commercial name	Bone injection needle beveled tip		
Manufacturer's reference	See table below		
Packaging	PE sachet in Tyvek pouch. To improve the safety of the sterile packaging, the cannula of the needles is protected with a PE tubing.		
Sterilization - Shelf life	Sterilization: Ethylene Oxide Shelf life: 3 years.		
Product description	Beveled tip with triple sharpening; graduated cannula; Luer-Lok connections suitable for biopsy purposes		
Composition	Cannula: steel AISI 304; stylet: steel AISI 302; handgrip: PP; Luer Lock cone and stylet knob: nickel-plated brass		
Quality controls	Raw materials: visual check, tensile and torsion tests. End product: visual inspection, label's verification, sterilization certificate check. Packaging: dimensional inspection, weld joint control (peeling test), leak test (methylene blue), control of weld thickness.		
Configuration	1 Bone Injection Needle		
M.D. classification	Ila according to CEE 93/42		
Intended use	Needle for acrylic bone cement injection into vertebral body.		
Latex content	The medical device and its packaging are latex-free		
Biocompatibility	The component of the device which comes in contact with the human body is made with stainless steel suitable for medical use (ISO 9626).		
Warnings and storage conditions	The device is sterile unless the packaging is open or damaged. The device is single-use and intended for use on a single patient. Do not re-sterilize or reuse. Re-sterilization may cause an infection risk to the patient.		
Disposal	All wastes should be considered surgical material and must be disposed as such at the end of surgery.		
Release of chemicals	N/A		
The state of the s			

Ref.	Description		
ASB0000	Beveled tip needle	G:09	L:15 cm
ASB0010	Beveled tip needle	G:11	L:12 cm
ASB0020	Beveled tip needle	G:13	L:12 cm
ASB0120	Beveled tip needle	G:15	L:12 cm

Release date: 18/12/2015



TECHNICAL SHEET

BONE INJECTION NEEDLES BEVELED TIPS — PLASTIC

Rev	Date	Description of changes	
00	30/09/08	New issue	
01	18/12/15	New format. Information updating	