

CORONARY

PERIPHERAL

MICROCATHETERS

NEUROVASCULAR

ANGIOGRAPHIC

ALL PRODUCTS




ASAHI Gaia® Series

Good choice for the need of high penetration while maintaining torqueability. Chronic Total Occlusion indicated. Guide wire with ACT ONE and low tip load to improve deflection and torqueability.




TIP LOAD



TIP RADIOPACITY 15 cm

COATING **HYDROPHILIC* WITH UNCOATED DISTAL TIP**
Hydrophilic 40 cm – uncoated distal tip

TIP SHAPE




Mini Pre-shape



ASAHI Gaia Second
TUBE CLIP



ASAHI Gaia First
TUBE CLIP



ACT ONE >



ASAHI Gaia Third
TUBE CLIP

Ordering Information

To order, please call us toll-free at 855-286-9473, or email us at customersupport@asahi-intecc.com. To see a complete listing of all available ASAHI products, go to our [Full Product Listing](#) page.

PRODUCT	CATALOG NO.	TIP SHAPE	LENGTH (CM)	DIAMETER (INCH(MM))	RADIOPAQUE SEGMENT (CM)	SP CC LE (C
ASAHI Gaia First	AHW14R007P	Pre- Shape	190	0.010/0.014 (0.26/0.36)	15	15
ASAHI Gaia First, 300 cm	AHW14R307P	Pre- Shape	300	0.010/0.014 (0.26/0.36)	15	15
ASAHI Gaia Second	AHW14R008P	Pre- Shape	190	0.011/0.014 (0.28/0.36)	15	15
ASAHI Gaia Second, 300 cm	AHW14R308P	Pre- Shape	300	0.011/0.014 (0.28/0.36)	15	15
ASAHI Gaia Third	AHW14R011P	Pre- Shape	190	0.012/0.014 (0.30/0.36)	15	15

PRODUCT	CATALOG NO.	TIP SHAPE	LENGTH (CM)	DIAMETER (INCH(MM))	RADIOPAQUE SEGMENT (CM)	SP CC LE (C
ASAHI Gaia Third, 300 cm	AHW14R311P	Pre- Shape	300	0.012/0.014 (0.30/0.36)	15	15

*Coated with SLIP-COAT® coating



Indications for Use

ASAHI PTCA Guide Wires are intended to facilitate the placement of balloon dilation catheters during percutaneous transluminal coronary angioplasty (PTCA) and percutaneous transluminal angioplasty (PTA). The ASAHI PTCA Guide Wires are not to be used in the neurovasculature.

[< Back to Coronary](#)

[BACK TO TOP ^](#)

ASAHI FULL PRODUCT LISTING

Select from our comprehensive portfolio of purpose-driven tools to master all your interventional cases.

FOUR CORE TECHNOLOGIES
FULL PRODUCT LISTING
CUSTOMER SERVICE
QUALITY POLICY

DOWNLOAD APP

PRODUCT LISTING >



CONTACT ASAHI INTECC

Toll-free: 855-286-9473
customersupport@asahi-intecc.com

CONTACT FORM >

Copyright © 2022 ASAHI INTECC USA, INC. All Rights Reserved | [PRIVACY POLICY](#) | [COOKIE POLICY](#)

Please read the terms of use for further details. Unless otherwise specified, all product names appearing in this website are trademarks or registered trademarks owned by or licensed to ASAHI. All other registered trademarks or trademarks are the property of their respective owners. No use of any ASAHI trademark, trade name, or trade dress in this website may be made without the prior written authorization of ASAHI, except to identify the product or services of the company.