

# PRE-ASSEMBLED EEG HEAD CAPS WITH STANDARD-THICKNESS ELECTRODES AND WIRES LOCATED EXTERNALLY

**REUSABLE**

**Packaging: box containing 1 piece.**

| CODE           | WIRES    | SIZE     | Number of electrodes   | DESCRIPTION         | COLOUR   |
|----------------|----------|----------|------------------------|---------------------|--|
| <b>CUEE60S</b> | external | 58-62 cm | 21, standard thickness | Extra large, adult  | Blue      |
| <b>CUEE56S</b> | external | 54-58 cm | 21, standard thickness | Large, adult        | Red       |
| <b>CUEE52S</b> | external | 50-54 cm | 21, standard thickness | Medium-small, adult | Yellow    |
| <b>CUEE48S</b> | external | 46-50 cm | 21, standard thickness | Large, child        | Green     |
| <b>CUEE44S</b> | external | 42-46 cm | 21, standard thickness | Medium, child       | L. Blue   |
| <b>CUEE41S</b> | external | 38-42 cm | 21, standard thickness | Small, child        | Pink      |
| <b>CUEE38S</b> | external | 34-38 cm | 11, standard thickness | Large, neonatal     | Brown    |
| <b>CUEE34S</b> | external | 32-34 cm | 11, standard thickness | Medium, neonatal    | Orange  |
| <b>CUEE32S</b> | external | 30-32 cm | 11, standard thickness | Small, neonatal     | White   |



**CUSTOMIZED CONNECTION CABLES ARE AVAILABLE UPON REQUEST.**

YOU CAN ORDER EEG HEAD CAPS WITH DIRECT CONNECTION TO THE EEG BOX.  
JUST REPLACE THE PRODUCT CODE'S LAST LETTER (**S**) WITH THE LETTER STANDING FOR THE  
DESIRED CONNECTION. SEE THE CHART BELOW.

|          |               |  |           |             |
|----------|---------------|--|-----------|-------------|
| <b>A</b> | INTERNATIONAL |  | <b>N</b>  | NEUROSOFT   |
| <b>M</b> | MICROMED      |  | <b>XX</b> | CUSTOM-MADE |