

TC-C3 Compact Trolley



Compact Trolley

TC-C3

The TC-C3 Compact Trolley is designed as a compact trolley to be used in conjunction with Optera CV-170.

Compact

The compact design means the TC-C3 is ideal for use in the outpatient or doctor's office where space is limited.

Easy to Move

The Compact Trolley is easy to move and position due to its light weight and ergonomically designed handles at the front of the top tray.

Easy to Clean

The TC-C3 is easy to clean due to limited dirt traps and no exposed fittings.

Monitor Support

The TC-C3 is offered with a monitor mount that supports monitors up to 26" with unrestricted swivel and tilt movement, allowing good viewing positions of the procedure image.

Large Top Tray

The top tray provides a clear working space at a comfortable height. The increased range of accessories enables the trolley to be utilized for a range of procedures.

Integrated Separation Transformer

The TC-C3 can be supplied with a separating transformer that is placed under the base, this means increased shelf space for equipment loading.



Accessories

TC-C3

MAJ-2146

Sliding Keyboard Tray



MAJ-1684

Scope Pole Kit



MAJ-2148

Rear Panel



MAJ-1686

Drawer Unit



MAJ-1685

IV Bag Hanger



MAJ-1683

LCD Monitor Mount



MAJ-1656

Footswitch Holder



MAJ-1639

CO₂ Cylinder Holder (140 mm)



MAJ-1665

Camera Head Holder



MAJ-2145

ENT Scope Hanger



MAJ-1658

Universal Storage Container



MAJ-110

Jar Mounting Kit TC



TC-C3

Standard sets include:

TC-C3 Compact Trolley plus items marked (•)

Accessories available as an optional extra (o)

			Region	Set 1	Set 2	Set 3	Set 4
MAJ	Part	Part	EUROPE				
Description	Description	Number	UK				
			ROW	K10030123	K10034675	K10030133	K10034676
				K10030124	K10034675	K10030134	K10034676
				K10030125	K10034675	K10030135	K10034676
MAJ-2146	Sliding Keyboard Tray	K10027573	o	•	•	•	•
MAJ-1684	Scope Pole Kit	K10027281	o	•	•		
MAJ-2148	Rear Panel	K10027574	o				
	Transformer			•		•	
MAJ-1686	Drawer Unit	K10027569	o				
MAJ-1685	IV Bag Hanger	K10027568	o				
MAJ-1683	LCD Monitor Mount	K10026888	o	•	•	•	
MAJ-1656	Footswitch Holder	K10021790	o				
MAJ-1639	CO ₂ Cylinder Holder (140 mm)	K10021041	o				
MAJ-1665	Camera Head Holder	K10028704	o				
MAJ-2145	ENT Scope Hanger	K10026497	o				
MAJ-1658	Universal Storage Container	K10021792	o				
MAJ-110	Jar Mounting Kit TC	K7503556	o				
MAJ-1654	Equipotential Terminal Strip	K10021352	o				
	Blister Pack Toolkit		o	•	•	•	•
	IEC Lead Set		o	•		•	

Specifications

Dimensions	1,110 mm (H) × 563.9 mm (D) × 489.5 mm (W)
	Height to top tray work surface – 1,100 mm
	Height with 26" LCD monitor in highest position – 1,945 mm
Weight	Maximum 82 kg (unladen) including fitted transformer

Voltage

The TC-C3 can be supplied with a separating transformer based on one of the following voltage options:

MAJ-2205	Transformer 110-120 V
MAJ-2206	Transformer 220-240 V

Trolley Construction

Columns	Extruded profiled aluminum
Shelves	Mild steel
Top Tray	Injection-moulded ABS

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG
 Postbox 10 49 08, 20034 Hamburg, Germany
 Wendenstrasse 14-18, 20097 Hamburg, Germany
 Phone: +49 40 23773-0, Fax: +49 40 233765
 www.olympus-europa.com