

## **.:** KEY FEATURE :.

Current output gas sensor that shows high signal stability at low cross interferences to anesthesia gases combined with a superior linearity over the entire range.

All characteristics are based on conditions at 25°C, 50% RH and 1013 hPa.

Measurement Range:	0 to 100 Vol.%
Expected Operating Life:	> 300,000 Vol.% h
Electrical Connector:	gold-plated slip rings on PCB
Initial Output Signal:	200 $\mu\text{A}\pm$ 40 $\mu\text{A}$ @ dry ambient air
Response Time t <sub>90</sub> :	< 20 s
Drift:	< 2 % Vol. O <sub>2</sub> /month @ air, averaged across 12 months
Operating Temperature:	10 to 40 °C
Pressure Range:	700 to 1250 hPa
Linearity Error:	$\leq$ 3 % @ 100% O <sub>2</sub> , applied for 5 min
Zero Offset:	less then 0.4% of oxygen equivalent @ 25 $^\circ\text{C}$ in 100% $N_2$ applied for 3 min
Repeatability:	$\pm$ 1 % Vol. O_2 @ 100 % O_2, applied for 5 $$ min
Influence of Humidity:	- 0.03 % rel. O <sub>2</sub> reading per % RH
Interferences:	according to DIN EN ISO 80601-2-55
Weight:	approximately 32 g
Material in Contact with Media:	PVDF, PPS, PTFE, stainless steel



## .: STORAGE CONDITIONS :.

**Temperature Range:** 

Humidity: Shelf Life: recommended: 5 to 30 °C maximum: - 20 to 50 °C (≤10h) up to 100 % RH, non-condensing < 3 months recommended

## **.:** RELATED PRODUCTS :.

Product O<sub>2</sub> - Sensor M-05 **Part-No.** 41 00 07 Housing Colour grey **Other Specifics** flat base

This data sheet is subject to change without prior notice. [M-05\_Rev04-2019\_0731.doc.doc]

IT Dr. Gambert GmbH :: Hinter dem Chor 21 :: 23966 Wismar :: Germany Phone: +49 (0)3841 220 050 :: Fax: +49 (0)3841 220 052 2 :: E-Mail: sales@itg-wismar.de

www.itg-wismar.de