



MEDJET

Sterile single use
foldable lens injectors



MEDJET | Sterile single use foldable lens injectors

Single use injectors for 1-piece foldable intraocular lenses

MEDJET MC^{1.6}	Type	Injector for hydrophilic IOLs
	Estimated incision size	1.6–1.7 mm
	Packaging	Cartridge + injector in blister; 1 disposable injection kit / box

MEDJET MB^{1.8}	Type	Injector for hydrophilic IOLs
	Estimated incision size	1.8–2.0 mm
	Packaging	Cartridge + injector in blister; 1 disposable injection kit / box

MEDJET MA^{2.2}	Type	Injector for hydrophilic IOLs
	Estimated incision size	2.2–2.3 mm
	Packaging	Cartridge + injector in blister; 1 disposable injection kit / box

MEDJET B1B^{2.2}	Type	Injector for hydrophilic IOLs
	Estimated incision size	2.2–2.3 mm
	Packaging	Cartridge + injector in peel pouch; 1 single use injector / box

MEDJET MX^{2.4} HB	Type	Injector for hydrophobic IOLs
	Estimated incision size	2.4–2.5 mm
	Packaging	Cartridge + injector in blister; 1 disposable injection kit / box

Technical data

Pusher	Silicon or TPE cushion for hydraulic injections
Sterilization	Ethylene Oxide

Cartridge

Winglets	with special click-lock system
Loading chamber	with guiding rails for an easy placement of the lens

EA-IMCMJEN 08 201504