

dispoSULIN insulin pen needle, sterile TD-42

Attachment No. 3 point 2 pp. 2.3 Revision I

2016.04.27

Page 1/1

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- For subcutaneous injection
- Dedicated for use with pre-filled pens (insulin pens)

FEATURES:

- Made of polypropylene, polyethylene and stainless steel
- The triple-point needle blade allows easy and atraumatic skin penetration
- Covering the needle with silicone lubricant reduces resistance during injection
- Thin walls of the needle ensure easy flow of the drug and the right flow rate
- The colored cap and the paper cover of the packaging make it easy to quickly recognize the size of the needle
- Available in various sizes, allowing you to adapt to individual needs
- Universal, compatible with all pen models
- Single use
- Non-pyrogenic
- Latex-free
- DEHP-free
- EO sterilized
- Shelf life: 5 years
- Intermediate packaging: 100 pcs



SIZE AND REFERENCE					
CATALOGUE NUMBER	SIZE	COLOUR	G	INTERMEDIATE PACKAGING	BULK PACKAGING
INP-033X12	0,33 mm x 12 mm	pink	29	100 pcs	50 x 100 pcs
INP-030X08	0,30 mm x 8 mm	yellow	30	100 pcs	50 x 100 pcs
INP-025X06	0,25 mm x 6 mm	violet	31	100 pcs	50 x 100 pcs
INP-025X08	0,25 mm x 8 mm	blue	31	100 pcs	50 x 100 pcs