PhysioGo 100A / 101A -

electrotherapy unit





ERGONOMICS

- two completely independent treatment channels
- · operation mode: program / manual
- preset treatment program and sequence base
- statistics
- encyclopedia with treatment methodology description
- · treatment programs chosen by name or medical field
- · concise list of favorite programs and sequences
- buzzer volume setting
- possibility of edition the names of user-defined programs and sequences
- self-test systematic control of the unit
- 7" color touchscreen
- · real time clock
- model PHG101A battery unit



ELECTROTHERAPY

- operation in CC (stabilized output current) and CV (stabilized output voltage) modes
- · electrodes test
- generation of unidirectional (unipolar) currents in interrupted mode
- complete galvanic isolation between channels in each mode

CURRENTS

- interferential (dynamic, isoplanar, one-channel (AMF))
- TENS (symmetric, asymmetric, alternating, burst)
- TENS to spastic paralysis therapy
- Kotz'/ Russian stimulation
- tonolysis
- diadynamic (MF, DF, CP, CP-ISO, LP)
- pulse rectangular, triangular
- Trabert pulse, Leduc' pulse
- unipolar sine surge
- galvanic
- microcurrents



ELECTRODIAGNOSTICS

- electrodiagnostics with graphic presentation of I/t curve
- automatic calculation of rheobase, chronaxie as well as threshold, coefficient and quotient of accommodation

MANUAL MODE

complete control over parameters for advanced users

TREATMENT PROGRAMS

- · 69 preset treatment programs
- 50 user-defined programs
- 69 favorite programs

TREATMENT SEQUENCES

- 38 preset treatment sequences
- 10 user-defined sequences



TECHNICAL SPECIFICATION

• controller – max constant current in patient's circuit (CC mode)

– galvanic	40 mA
 diadynamic, pulse 	60 mA
 interferential, Kotz, unipolar sine surge 	100 mA
- TENS	140 mA
tonololysis	100 mA
microcurrents	1000 μΑ
 max voltage in patient's circuit (CV mode) 	140 V

mains supply
 230 V, 50/60 Hz, max 75W, 90VA

• dimensions 34 x 28 x 11-16 cm

• weight 6 kg



STANDARD ACCESSORIES

- mains cable
- patient's cable (2 pcs)
- electrotherapy electrodes 6 x 6 cm (4 pcs); 7,5 x 9 cm (2 pcs)
- viscose electrode covers 6 x 6 cm (8 pcs); 7,5 x 9 cm (4 pcs)
- velcro fixing belts 40 x 10 cm (2 pcs); 100 x 10 cm (2 pcs)
- spare fuses WTA-T1 A/250 V (2 pcs)
- touchscreen stylus pen (1 pcs)
- touchscreen cleaning cloth (1 pcs)

OPTIONAL ACCESSORIES

- disposable self-adhesive electrodes 3 x 4 cm; 5,5 x 5,5 cm; 5,5 x 12,5 cm
- point electrodes 6 mm, 10 mm, 15 mm, 20 mm
- · crocodile clip
- sand-filled bag 21 x 14 cm, 21 x 28 cm
- combination therapy cable
- patient's stop switch
- bag for the unit and accessories
- Versa trolley
- battery

