



elastopor IV
**IV cannula dressing, non-woven,
self-adhesive, sterile**

TD-09

att.no.3 pt.2 subpt.2.3

revision 2

date: 2022-01-05

amendment of: 2023-02-01

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- stabilization, protection and fixation of intravenous cannulas

FEATURES:

- self-adhesive dressing with a notch for the vertical port of the intravenous cannula
- made of hydrophobic non-woven
- centrally situated absorbent pad with polyethylene mesh protects the non-woven from sticking to the injection site
- non-woven lining should be placed under cannula's wings, which prevents it from causing damage to the skin
- hypoallergenic
- rounded edges prevent accidental removal of the dressing
- coated with a hypoallergenic acrylic adhesive without zinc oxide, applied to the non-woven fabric by a continuous method
- silicone paper backing facilitates precise, painless and efficient application
- graphical instructions for use on individual and intermediate packaging
- EO sterilized
- shelf life: 5 years



SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	SALES PACKAGING BOX	TRANSPORT PACKAGING
811013E	6 cm x 8 cm	1 pc.	100 pcs	20 x 100 pcs
811012E	5,1 cm x 7,6 cm	1 pc.	100 pcs	20 x 100 pcs
811023E	6 cm x 8 cm	1 pc.	50 pcs	40 x 50 pcs
811022E	5,1 cm x 7,6 cm	1 pc.	50 pcs	40 x 50 pcs