



elastoDERM F-IV
cannula dressing, foil, with frame,
with U-incision, self-adhesive, sterile

att.no.3 pt.2 subpt.2.3

revision 1

date: 2023-10-08

amendment of: N/A

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- fixing and protecting intravenous cannulas and protect and observe the insertion site

FEATURES:

- Self-adhesive polyurethane film dressing having a cutout for the vertical port of a peripheral intravenous cannula
- rounded edges prevent the dressing from accidentally coming off
- transparent - allows observation of the application site without removing the dressing
- vapor - permeable, ensuring proper gas exchange between the dressing and the skin, which is protected from excessive moisture accumulation and maceration
- provides a barrier to microorganisms and bacteria minimizing the risk of infection
- thin and elastic - adheres well, adapts to the shape of the body
- does not absorb X-rays
- acrylic adhesive applied continuously over the entire surface of the adhesive
- has a protective layer of siliconized paper to facilitate precise, painless and effective application
- wings reinforced from the inside with a layer of nonwoven fabric
- presence of two fixing non-woven strips
- has a "frame" type application system to facilitate precise and effective application of the dressing
- a descriptive tape allowing you to note, for example, the date the dressing was applied
- skin contact up to 7 days
- waterproof
- hypoallergenic
- single-use
- sterile
- shelf life: 3 years



SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	SALES PACKAGING BOX	TRANSPORT PACKAGING
812009	9 cm x 7,2 cm	1 pc.	100 pcs	24 x 100 pcs