

elastopor STERIL D non-woven dressing with absorbent pad, incision and O-hole, self-adhesive, sterile

TD-09

att.no.3 pt.2 subpt.2.3 revision 2

date: 2022-09-18

amendment of: 2023-02-01

MANUFACTURER:

ZARYS International Group sp. z o.o. sp.k.

ADDRESS:

ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- protection of the site of drain or catheter application
- dressing wounds with small and medium exudate level

FEATURES:

- self-adhesive dressing with absorbent pad
- dressing with incision and central O-hole
- incision and O-hole enabling application around a fixed drain
- made of hydrophobic non-woven material
- micropores in the non-woven material structure provide an adequate rate of moisture vapour permeability so correct gas exchange between the dressing and skin is guaranteed
- centrally situated absorbent pad absorbs small amounts of exudates
- absorbent pad coated with polyethylene mesh does not stick to the wound as well as facilitates removal of the dressing
- hypoallergenic
- elastic comfortable for the patient, adjusts to the body shape and guarantees freedom of movement
- easy to attach to the skin
- rounded edges prevent accidental removal of the dressing
- covered with hypoallergenic acrylic glue
- silicone paper backing facilitates precise, painless and efficient application
- EO sterilized
- shelf life: 5 years



SIZE AND REFERENCE					
CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	SALES PACKAGING BOX	TRANSPORT PACKAGING	
801031E	9 cm x 10 cm	1 pc.	30 pcs	24 x 30 pcs	
801030E	12 cm x 14 cm	1 pc.	25 pcs	16 x 25 pcs	

Marketing material intended for the health care professionals. Prepared by: W.W. – Product Manager ZARYS International Group sp. z o.o. sp.k.	page 1/1
No part of this document may be changed or copied without the permission of ZARYS International Group sp.z o.o. sp.k.	rev.2023.01