

# SRI TRANG GLOVES (THAILAND) PUBLIC COMPANY LIMITED

10 Soi 10 Phetkasem Road, Hat Yai, Hat Yai, Songkhla 90110 Thailand  
Tel. (66) 074-344-663 Fax. (66) 074-344-677

SCT.PD.SP.05.035-010622 R.05

PRODUCT SPECIFICATION SHEET : LOOOGF-B-EU-M-NS

DAR-SCT-2022-0614

PAGE 1 of 1

## GLOVE SPECIFICATION

**Product Specification Group Code :** LOOOGF-B-EU-M-NS  
**Product Code :** DLOFBOG  
**Product name :** Latex Examination Gloves, Powder Free, Chlorinated, Non Sterile  
**Grade :** EU Spec, Medical Grade

### GENERAL INFORMATION

**Type :** Single-use disposable examination or protective glove  
**Marking :** Gloves are not marked to designated size  
**Shape :** Ambidextrous  
**Material :** Natural Rubber  
**Warning:** This product contains natural rubber latex which may cause allergic reactions  
**Cuff :** Rolled cuff  
**Surface :** Textured at Finger Tip

<b>Thickness :</b> (single wall according to ASTM D3578, ASTM D3767)	<b>Finger :</b> Min. 0.08mm	<b>Palm :</b> Min. 0.08mm
---	-----------------------------	---------------------------

<b>Residual powder :</b> (according to EN 455-3, EN ISO 21171)	Max. 2.0mg/glove
--	------------------

<b>Force at break :</b> (according to EN 455-2)	Before Aging : Med. $\geq$ 6.0N After Aging : Med. $\geq$ 6.0N
--	---

<b>Protein content :</b> (according to EN 455-3)	Max. 50 $\mu$ g/g
---	-------------------

<b>Tensile Strength :</b> (according to ASTM D3578, ASTM D412, ASTM D573)	Before Aging : Min. 18Mpa After Aging : Min. 14Mpa
---	---

<b>Length :</b> (according to EN 455-2)	XS-XL : Med. $\geq$ 240mm
--	---------------------------


<b>% Elongation :</b> (according to ASTM D3578, ASTM D412, ASTM D573)	Before Aging : Min. 650% After Aging : Min. 500%
---	---

<b>Available Sizes &amp; Palm width (mm) :</b> XS (5-6) : $\leq$ 80, S (6-7) : 80+/-10, M (7-8) : 95 +/-10, L (8-9) : 110+/-10, XL (9-10) : $\geq$ 110 (according to EN 455-2)
---

**Storage:** 5 years storage life (date of expiration see on the box). Keep area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper-ions discolor the glove. Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead to accelerated aging and should be avoided.

### STANDARDS, GUIDELINES AND QUALITY CERTIFICATION

(Every mentioned standard is used in the latest edition.)

<b>Certification :</b>	ISO 9001, ISO 13485															
<b>Conformity to guidelines :</b>	Medical Device Regulation (EU)2017/745 : Class I under Rule1 and Rule 5 according to Annex VIII (Subject to EU Declaration of Conformity to brand owner) Personal Protective Equipment Regulation (EU) 2016/425 : Category I or Category III (Subject to EU Type-Examination Certificate)															
<b>Conformity to standards :</b>	<table><tr><td><u>Medical</u></td><td><u>PPE</u></td><td><u>PPE</u></td></tr><tr><td>EN 455</td><td>CAT I</td><td>CAT III</td></tr><tr><td>ASTM D3578</td><td>EN ISO 21420</td><td>EN ISO 21420</td></tr><tr><td></td><td></td><td>EN ISO 374</td></tr><tr><td></td><td></td><td>EN 16523-1</td></tr></table>	<u>Medical</u>	<u>PPE</u>	<u>PPE</u>	EN 455	CAT I	CAT III	ASTM D3578	EN ISO 21420	EN ISO 21420			EN ISO 374			EN 16523-1
<u>Medical</u>	<u>PPE</u>	<u>PPE</u>														
EN 455	CAT I	CAT III														
ASTM D3578	EN ISO 21420	EN ISO 21420														
		EN ISO 374														
		EN 16523-1														
<b>Conformity to food regulation :</b>	1935/2004 															

### WARNING AND REPORTING SYSTEM

**Potentially allergenic ingredients :** Dithiocarbamate type.  
(For further information, the actual version of the ingredients list is given on request.)

Please report any incidents with the products to **STGT\_TA@sritranggroup.com** and **stgt\_tp@sritranggroup.com** without undue delay. Our Product Safety Officers will review all information provided and might contact the original informant for further clarification.

## SRI TRANG GLOVES (THAILAND) PUBLIC COMPANY LIMITED

10 Soi 10 Phetkasem Road, Hat Yai, Hat Yai, Songkhla 90110 Thailand  
Tel. (66) 074-344-663 Fax. (66) 074-344-677

SCT.PD.SP.05.035-010622 R.05

PRODUCT SPECIFICATION SHEET : LOOOGF-B-EU-M-NS

DAR-SCT-2022-0614

Appendix of specification

### APPENDIX : 2

**Product Specification Code :** LOOOGF-B-WHT-EU-MC-NS  
**Product Code :** DLOFBOG  
**Material :** Natural Rubber  
**Treatment and Feature :** Powder Free, Online Chlorinated

### SAMPLING, INSPECTION AND FINAL RELEASE INFORMATION

Product Specification Code	Glove color	Characteristics			Test-Method
		Barrier Defects*	Major Visual Defects*	Minor Defects*	
LOOOGF-B-WHT-EU-MC-NS	WHITE (PALE YELLOW)	AQL ≤ 1.5	AQL ≤ 1.5	AQL ≤ 4.0	EN 455-1

#### **Barrier Defects :**

Inspection level G I (ISO 2859-1)

Holes that results in the malfunction of glove as a protection barrier either physically or not physically affect the appearance.

#### **Major Visual Defects :**

Inspection level G I (ISO 2859-1)

Defined defect that disrupt the glove's function or defect that cause an uncomfortable feeling to the user.

#### **Minor Defects :**

Inspection level G I (ISO 2859-1)

Visual defects aggregated are non-conformities of lower degree of concern, which do not prevent correct use of gloves. Minor cosmetic defect which does not affect on the glove's performance

Remark :

\* Defect Characteristics following In-house practice.

<b>Product description / Details of packaging for Change Control No.:</b> OTC-STA-0821-1675			
<b>1. Brand Name:</b> SRITRANG		<b>Customer :</b> IMUNOTEHNOMED(MOLDOVA)	
<b>Brand Name (INT):</b> STGEUVVH-10X100-LOFBOG-EUC-WHT-NS-BX			
<b>2. Product Type:</b> LATEX POWDER FREE WHITE TEXTURED AT FINER TIPS NON-STERILE			
<b>3. KMAT :</b> DLOFBOG		<b>510K No.:</b> -	
<b>Specification:</b> LOOOGF-B-WHT-EU-MC-NS		<b>CA No.:</b> -	
<b>4. Material number (FG):</b>			
Size XS : 14019467		Size L : 14019470	
Size S : 14019468		Size XL : 14019471	
Size M : 14019469		Size 2L : -	
<b>5. Packaging</b>			
<b>Packaging type :</b> INNER BOX			
<b>Dimension(mm):</b> PKGPF02    W: 122    L: 240    H: 65			
SIZE	COLOR	MATERIAL NO.(PCK)	VERSION
FS	BLACK,369C	05012002301304443 (20014203)	C00
XS			
S			
M			
L			
XL			
2L			
<b>Packaging type :</b> OUTER BOX			
<b>Dimension(mm):</b> PKGPF02    W: 249    L: 340    H: 250			
SIZE	COLOR	MATERIAL NO.(PCK)	VERSION
FS	BLACK,369U	05011002301303744 (20014204)	C00
XS			
S			
M			
L			
XL			
2L			
<b>Packaging type :</b>			
<b>Dimension(mm):</b> W:    L:    H:			
SIZE	COLOR	MATERIAL NO.(PCK)	VERSION
FS			
XS			
S			
M			
L			
XL			
2L			
<b>Packaging type :</b>			
<b>Dimension(mm):</b> W:    L:    H:			
SIZE	COLOR	MATERIAL NO.(PCK)	VERSION
FS			
XS			
S			
M			
L			
XL			
2L			
<b>Packaging type :</b>			
<b>Dimension(mm):</b> W:    L:    H:			
SIZE	COLOR	MATERIAL NO.(PCK)	VERSION
FS			
XS			
S			
M			
L			
XL			
2L			

<b>Material:</b> DUPLEX	<b>FSC:</b> NO
<b>Flap:</b> LOCKING	<b>Coating:</b> VARNISHING
<b>Opening:</b> OTHER AND NON STD-SINGLE LINE	
<b>Marking Size by:</b> TICK SIZE BY BLACK	
<b>MFG+Lot No.:</b> YYYY-MM(MFG) XXXX(INT) CCCC(CASE) DLOFBOG(SPACIAL)	
<b>MFG:</b> -	
<b>EXP :</b> YYYY-MM(MFG+3YRS)	
<b>Pre-print:</b> NO	
<b>Remark:</b> Share artwork with : TGEUVVH-10X100-LOFBOG-EUC-WHT-NS-BX Inkjet "DLOFBOG" behind (CASE NO.)	
<b>Material:</b> B,BROWN	<b>FSC:</b> NO
<b>Flap:</b> -	<b>Coating:</b> -
<b>Opening:</b> -	
<b>Marking Size by:</b> STAMP TEXT BY BLACK	
<b>MFG+Lot No.:</b> YYYY-MM(MFG) XXXX(INT) CCCC(CASE) REF:DLOFBOG(SPACIAL)	
<b>MFG:</b> -	
<b>EXP :</b> YYYY-MM(MFG+3YRS)	
<b>Pre-print:</b> NO	
<b>Remark:</b> Share artwork with : TGEUVVH-10X100-LOFBOG-EUC-WHT-NS-BX Inkjet REF:DLOFBOG short side 2	
<b>Material:</b>	<b>FSC:</b>
<b>Flap:</b>	<b>Coating:</b>
<b>Opening:</b>	
<b>Marking Size by:</b>	
<b>MFG+Lot No.:</b>	
<b>MFG:</b>	
<b>EXP :</b>	
<b>Pre-print:</b>	
<b>Remark:</b>	
<b>Material:</b>	<b>FSC:</b>
<b>Flap:</b>	<b>Coating:</b>
<b>Opening:</b>	
<b>Marking Size by:</b>	
<b>MFG+Lot No.:</b>	
<b>MFG:</b>	
<b>EXP :</b>	
<b>Pre-print:</b>	
<b>Remark:</b>	



St.B.



## 6. Packing/Contain instructions.

Dispenser Approx. 100(XL:90) pcs by weight/dispenser box.

☐ Use ☒ Not Use Plastic under the perforation of dispenser boxes

Carton 10 Dispenser boxes./ Carton.

☐ Use ☒ Not Use Plastic bag for dispenser boxes

Packing Method LAYER

☐ Use ☒ Not Use Poly bag for carton Color:

## 7. Loading quantity

Container 20ft : 1392 Cartons / Container 40ft : 2838 Cartons / Container 40ftHC : 3220 Cartons / Truckt : Cartons

## 8. Change description



For Dispenser box :

Add customer IMUNOTEHNOMED(MOLDOVA) under brand :  
 STGEUVVH-10X100-LOFBOG-EUC-WHT-NS-BX  
 Share artwork with brand STGEUVVH-10X100-LOFBOG-EUC-WHT-NS-BX  
 (NEXTECH(TURKEY) CUSTOMER)


For Transport carton :

Add customer IMUNOTEHNOMED(MOLDOVA) under brand :  
 STGEUVVH-10X100-LOFBOG-EUC-WHT-NS-BX  
 Share artwork with brand STGEUVVH-10X100-LOFBOG-EUC-WHT-NS-BX  
 (NEXTECH(TURKEY) CUSTOMER)

## 9. Packing format and Inkjet format



INKJET BY STGT INTERNAL DEPT.



St.B.