

STAMMBERGER Circular Cutting Punches, with cleaning connector

Special Features:


- For enlarging openings in the sphenoid frontal wall
- Circular cutting punch mechanism allows cutting in a full circle of 360° without rotating the instrument as required with a conventional punch forceps
- No interference with other instruments simultaneously used in the nose (e. g. endoscope, suction tube)
- Available in 2 sizes: diameter 3.5 and 4.5 mm, punch head 4-fold LASER-welded
- With cleaning connector
- Multipurpose use:
 - In addition to abrading procedures applied to the sphenoid frontal wall, bony ethmoid septa, pieces of nasal concha and other thin bony bridges can also be cut away
 - Suitable for treatment of choanal atresia



651055 STAMMBERGER **Punch**, circular cutting, for sphenoid, ethmoid and choanal atresia, diameter 3.5 mm, with cleaning connector, working length 18 cm



651050 **Same**, diameter 4.5 mm



651060 STAMMBERGER **Punch**, circular cutting, 65° upturned, for frontal sinus recess, diameter 3.5 mm, with cleaning connector, working length 17 cm

651065 **Same**, diameter 4.5 mm