


Foreign Body Holding Forceps

Used to retrieve foreign bodies from gastrointestinal tract

Shark Tooth Forceps can be used with ERCP scope to retrieve Biliary-Stent



Art. No.	Productname	Required Channel Ø (mm)	Jaws Specification	Working Length (cm)
SSTF-170	Shark Tooth Forceps	2.8	Rat Tooth	170
SSTF-230	Shark Tooth Forceps	2.8	Rat Tooth	230
SATF-170	Alligator Forceps	2.8	Alligator	170
SATF-230	Alligator Forceps	2.8	Alligator	230
SSTF-170-R	Shark Tooth Forceps (Rotatable)	2.8	Rat Tooth	170
SSTF-230-R	Shark Tooth Forceps (Rotatable)	2.8	Rat Tooth	230
SATF-170-R	Alligator Forceps (Rotatable)	2.8	Alligator	170
SATF-230-R	Alligator Forceps (Rotatable)	2.8	Alligator	230
SSAF-170	Shark Tooth Alligator Combo	2.8	Combo	170
SSAF-230	Shark Tooth Alligator Combo	2.8	Combo	230
SSAF-170-R	Shark Tooth Alligator Combo (Rotatable)	2.8	Combo	170
SSAF-230-R	Shark Tooth Alligator Combo (Rotatable)	2.8	Combo	230

Class IIa device – non sterile – reusable

Memory Basket (Nitinol)

Used to retrieve foreign bodies/stones from gastrointestinal tract through the scope

Not for Lithotripsy

Basket Type: Oval, Hexagonal, Tipless and Helical
Retains shape for longer time
Light handle with flush port for contrast liquid



Art. No.	Basket Ø (mm)	Working Length (cm)	Working Channel Ø (mm)
SMMB-15-T	15	230	2.8
SMMB-20-T	20	230	2.8
SMMB-25-T	25	230	2.8
SMMB-30-T	30	230	2.8
SMMB-35-T	35	230	2.8
SMMB-40-T	40	230	2.8
SMMB-45-T	45	230	2.8

T = Tipless

For other shapes please add to Art. No. -O = Oval, -HX = Hexagonal and -H = Helical