


RADIFOCUS™ *Introducer II*

TRANSRADIAL KIT


Introducer Sheath



Introducers are intended to be inserted percutaneously into a vessel to facilitate the entire interventional procedure. Being a pioneer and a leader in vascular access, Terumo's goal is to offer complete solutions for transradial interventions with devices that minimize patient stress and optimize outcomes.

Product Characteristics

- Dilator is tapered to each guidewire to minimize vessel trauma
- Terumo cross-cut valve maintains haemostasis to avoid bleeding and air aspiration
- Thin radiopaque sheath with anti-kinking sleeve for easy catheter handling
- Snap-on/click-off dilator to prevent dilator back-out during insertion and allows one-hand unlocking
- Micro puncture metallic entry needle with short bevel (22G / 21G / 20G) equals to mini guidewire to minimize puncture site complications



- A Mini guidewire
B Dilator
C Sheath
D Fr size indication
E Side tube
F 3 way stop cock

General specifications

Mini Guidewire - Type	Spring Straight
Mini Guidewire - Length	45 cm
Entry Needle - Type	Metal
Entry Needle - Length	35 mm
Packaging	Tray
R Kit content	Sheath, dilator, spring mini guidewire and metallic entry needle

Item specifications

Sheath Size	Sheath Length	Dilator Tip - Inner Diameter	Entry Needle - Outer Diameter	Guidewire Compatibility - Diameter	Code
4 Fr	7 cm	0.018 in min. 0.49 mm	22 G 0.7 mm	0.018 in max. 0.47 mm	RT-R40A07PQ
4 Fr	7 cm	0.025 in min. 0.65 mm	20 G 0.9 mm	0.025 in max. 0.64 mm	RT-R40G07PQ
4 Fr	10 cm	0.018 in min. 0.49 mm	22 G 0.7 mm	0.018 in max. 0.47 mm	RT-R40A10PQ
4 Fr	10 cm	0.021 in min. 0.56 mm	21 G 0.8 mm	0.021 in max. 0.54 mm	RT-R40D10PQ
4 Fr	10 cm	0.025 in min. 0.65 mm	20 G 0.9 mm	0.025 in max. 0.64 mm	RT-R40G10PQ
5 Fr	7 cm	0.018 in min. 0.49 mm	22 G 0.7 mm	0.018 in max. 0.47 mm	RT-R50A07PQ
5 Fr	7 cm	0.025 in min. 0.65 mm	20 G 0.9 mm	0.025 in max. 0.64 mm	RT-R50G07PQ
5 Fr	10 cm	0.018 in min. 0.49 mm	22 G 0.7 mm	0.018 in max. 0.47 mm	RT-R50A10PQ
5 Fr	10 cm	0.021 in min. 0.56 mm	21 G 0.8 mm	0.021 in max. 0.54 mm	RT-R50D10PQ
5 Fr	10 cm	0.025 in min. 0.65 mm	20 G 0.9 mm	0.025 in max. 0.64 mm	RT-R50G10PQ
6 Fr	7 cm	0.018 in min. 0.49 mm	22 G 0.7 mm	0.018 in max. 0.47 mm	RT-R60A07PQ
6 Fr	7 cm	0.025 in min. 0.65 mm	20 G 0.9 mm	0.025 in max. 0.64 mm	RT-R60G07PQ
6 Fr	10 cm	0.018 in min. 0.49 mm	22 G 0.7 mm	0.018 in max. 0.47 mm	RT-R60A10PQ
6 Fr	10 cm	0.021 in min. 0.56 mm	21 G 0.8 mm	0.021 in max. 0.54 mm	RT-R60D10PQ
6 Fr	10 cm	0.025 in min. 0.65 mm	20 G 0.9 mm	0.025 in max. 0.64 mm	RT-R60G10PQ
7 Fr	10 cm	0.021 in min. 0.56 mm	21 G 0.8 mm	0.021 in max. 0.54 mm	RT-R70D10PQ

Unit per box: 5 pcs.

Please quote above item reference codes when placing an order