









EKG-10005B	K019 (BA-901D)	10 leads; integrated; AHA, DB15M-Din3.0 	 <p><i>Nihon Kohden ECG-9620, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9130, ECG-9132, DongJiang ECG-11A, ECG-11B, ECG-11D, ECG-32A, (II), ECG-120 For ECG-1150/1250/1350/2350</i></p>
EKG-10005C	K081 (BJ-902D)	10 leads; integrated; IEC, DB15M-Banana4.0 	
EKG-10005D		10 leads, integrated; AHA, DB15M-Banana4.0 	
EKG-10005E		10 leads, integrated; IEC, DB15M-Clip end 	
EKG-10005F		10 leads, integrated; AHA, DB15M-Clip end 	
EKG-10005G	K082 BJ-903D	10 leads; integrated; IEC, DB15M-snap 	
EKG-10005H		10 leads; integrated; AHA, DB15M-snap 	
EKG-10005I	BJ-904D; K098	5 leads integrated; IEC, din3.0	<i>for ECG-1950</i>
EKG-10005J		5 leads integrated; AHA, din3.0	
EKG-10005K		5 leads integrated; IEC, banana 