JRON .	27610001	calcuson ultrasonic generator set, for use with ultrasound lithotripsy probes and transducers, power supply 100-120 VAC/200-240 VAC, 50/60 Hz including: Mains Cord One-Pedal Footswitch Transducer Connecting Cable (transducer/generator) Protection Sleeve, for storage and sterilization of probes Cleaning Rod, for probes
	20300051	Suction Bottle, 0.5 I, sterilizable
	20300039	Bottle Cap, for Suction Bottle 20300051 LOT 11
	UP011	Bottle Holder, for Suction Bottle 20300051

Further accessories can be found in the Urology catalog.