

(1) Multi Functions

- ClearCoajet is a combination of Injection catheter and Hemostat.
- Marking, Injection, Coagulation are all possible with this device for EMR, ESD procedure.

(2) Back-cut Needle design

• It has advantages to Insertion&injection saline because the shape of Needle is sharp and I.D(Inner Diameter) is large.




Injection without



SPECIFICATION

Model	Needle O.D	Needle Thickness	Needle Length	Sheath 0.D	Working Length	Number of Use	Validity of Sterilzation
FM-CJ0002(1800) FM-CJ0002(2300)	23G	0.1	4mm	¢2.4mm	1,800mm 2,300mm	Single Use Only	3 years

ClearHemograsper



(1) Rotation with the handle

- Enhanced the Rotation, and the users can perform Hemostasis for bleeding at demanding location more easily.
- Jaw and Sheath are rotated separately.

(2) Installed irrigation port

• Provide convenience to users that use the scope without Water Jet function.

Irrigation Function Syringe or water pump irrigation port

SPECIFICATION

Model	Opening width	Sheath O.D	Type of Rotation	Working Length	Number of Use	Validity of Sterilzation
FM-HG0001(1650) FM-HG0001(1820) FM-HG0001(2000)	5mm	¢2.3mm	Rotational	1,650mm 1,820mm 2,000mm	Single Use Only	3 years

Imagine What a Woderful difference You have 7