

## Clear-Therm™ range

### Clear-Therm Midi

A low volume option for minimising dead space in anaesthesia.

Make an enquiry



| Box Qty.                   | 100                       | 50                        | 75                                                                  |
|----------------------------|---------------------------|---------------------------|---------------------------------------------------------------------|
| Luer lock port             | ✓                         | ✓                         | ✓                                                                   |
| Moisture loss              | 13.3mg H <sub>2</sub> O/L | 13.3mg H <sub>2</sub> O/L | 13.3mg H <sub>2</sub> O/L                                           |
| Calculated moisture return | 25.9mg H <sub>2</sub> O/L | 25.9mg H <sub>2</sub> O/L | 25.9mg H <sub>2</sub> O/L                                           |
| Filtration efficiency      | >99.9%                    | >99.9%                    | >99.9%                                                              |
| Resistance at 30L/min      | 0.9cm H <sub>2</sub> O    | 1.2cm H <sub>2</sub> O    | 1.2cm H <sub>2</sub> O                                              |
| Resistance at 60L/min      | 2.1 cm H <sub>2</sub> O   | 3.1 cm H <sub>2</sub> O   | 3.2cm H <sub>2</sub> O                                              |
| Compressible volume        | 34ml                      | 34ml*                     | 34ml*                                                               |
| Weight                     | 19g                       | 26g                       | 35g                                                                 |
| Connectors                 | 22F-22M/15F               | 22F-22M/15F               | 22F-22M/15F                                                         |
| Minimum tidal volume       | >100ml                    | >100ml                    | >100ml                                                              |
| Accessories                |                           | Fixed elbow**<br>1997000  | Superset™ catheter mount<br>with fixed elbow**<br>3504000 + 1997000 |

### Clear-Therm Mini and Micro

For paediatrics and neonates