

## Memory Basket (Nitinol)

Used to retrieve foreign bodies/stones from gastrointestinal tract through the scope

### Not for Lithotripsy

Basket Type: Oval, Hexagonal, Tipless and Helical  
Retains shape for longer time  
Light handle with flush port for contrast liquid



Art. No.	Basket Ø (mm)	Working Length (cm)	Working Channel Ø (mm)
<b>SMMB-15-T</b>	15	230	2.8
<b>SMMB-20-T</b>	20	230	2.8
<b>SMMB-25-T</b>	25	230	2.8
<b>SMMB-30-T</b>	30	230	2.8
<b>SMMB-35-T</b>	35	230	2.8
<b>SMMB-40-T</b>	40	230	2.8
<b>SMMB-45-T</b>	45	230	2.8

T = Tipless

For other shapes please add to Art. No. -O = Oval, -HX = Hexagonal and -H = Helical


Class IIa device – sterile – for single use only

## Memory Basket (Wire guided)

Used to retrieve foreign bodies/stones from gastrointestinal tract through the scope

### Not for Lithotripsy

Handle light weight with flush port for Contrast-Drug-Flush  
Retains shape for longer time  
Centrally located wire guide entrance  
Wire guide compatibility 0.035"



Art. No.	Basket Ø (mm)	Working Length (cm)	Working Channel Ø (mm)
<b>SMMB-WG-30</b>	30	230	2.8
<b>SMMB-WG-35</b>	35	230	2.8
<b>SMMB-WG-40</b>	40	230	2.8

Class IIa device – sterile – for single use only