



Raucocel®

Epistaxis packings



Epistaxis packings 100



Epistaxis packings 55

Fields of application

- Epistaxis packing 100: for haemostasis in the posterior nasal area and for postoperative use, e.g. after conchotomy
- Epistaxis packing 55: for haemostasis in the medial nasal area

Properties

- physiological blood coagulation based on the principle of autocoagulation
- time-saving and atraumatic application
- can be used to deliver solutions containing medication

Product composition

polyvinyl alcohol, drawstring made from surgical silk

To note

- thoroughly moisten packing prior to removal

Advice

- We also offer the following related accessories:
 - ◆ instrument wipe
 - ◆ headband

Raucocel® Epistaxis Packing, 100 mm

individually sealed and sterile, with thread

Length (mm)	Dimensions (mm)	REF	Shipping units
100	3.8 x 12.0	10143	10/240

Raucocel® Epistaxis Packing, 55 mm

individually sealed and sterile, with thread

Length (mm)	Dimensions (mm)	REF	Shipping units
55	3.8 x 12.0	10144	20/240

Raucocel® Accessories

Length (mm)	Dimensions (mm)	REF	Shipping units
Instrument wipe, individually sealed and sterile			
76	76 x 2	10313	20/480
headband, individually sealed			
175	40 x 5	11971	10/240

Not available in all countries