


Foreign Body Holding Forceps

Used to retrieve foreign bodies from gastrointestinal tract

Shark Tooth Forceps can be used with ERCP scope to retrieve Biliary-Stent



Art. No.	Productname	Required Channel Ø (mm)	Jaws Specification	Working Length (cm)
SSTF-170	Shark Tooth Forceps	2.8	Rat Tooth	170
SSTF-230	Shark Tooth Forceps	2.8	Rat Tooth	230
SATF-170	Alligator Forceps	2.8	Alligator	170
SATF-230	Alligator Forceps	2.8	Alligator	230
SSTF-170-R	Shark Tooth Forceps (Rotatable)	2.8	Rat Tooth	170
SSTF-230-R	Shark Tooth Forceps (Rotatable)	2.8	Rat Tooth	230
SATF-170-R	Alligator Forceps (Rotatable)	2.8	Alligator	170
SATF-230-R	Alligator Forceps (Rotatable)	2.8	Alligator	230
SSAF-170	Shark Tooth Alligator Combo	2.8	Combo	170
SSAF-230	Shark Tooth Alligator Combo	2.8	Combo	230
SSAF-170-R	Shark Tooth Alligator Combo (Rotatable)	2.8	Combo	170
SSAF-230-R	Shark Tooth Alligator Combo (Rotatable)	2.8	Combo	230

Class IIa device – non sterile – reusable

Foreign Body Holding Forceps

2 Prong, 3 Prong and 4 Prong

Used to retrieve foreign bodies from gastrointestinal tract



Art. No.	Prongs	Working Length (cm)	Working Channel Ø (mm)
SFBH-223	2	170	2.8
SFBH-218	2	170	2.0
SFBH-323	3	170	2.8
SFBH-318	3	170	2.0
SFBH-423	4	170	2.8
SFBH-223-C	2	230	2.8
SFBH-218-C	2	230	2.0
SFBH-323-C	3	230	2.8
SFBH-318-C	3	230	2.0
SFBH-423-C	4	230	2.8

Class IIa device – reusable