

aurolab aqueous drainage implant



MODEL
aadi 350

aadi, a non-valved glaucoma shunt indicated for glaucoma patients for not responding maximal medical therapy, failed trabeculectomy, neovascular, congenital and uveitic glaucoma

CE
0197

- *Made of permanent implantable grade silicone, a proven material for patient safety*
- *The silicone plate that measures 13mm convex radius with surface area of 350mm² and with 35mm straight silicone tube takes care of curvature of the globe and easy insertion*
- *The two fixation holes in the silicone plate can be used to suture the plate to the sclera using 8-0 or 9-0 aurodon sutures*



Key Features of aadi

Larger surface Area

Larger surface area of the silicone plate provides aqueous filtration and pressure relief by controlling long term intraocular pressure with fewer post operative complications and medications

Single Quadrant Insertion

Innovative design for simple and fast surgical technique which allows surgeons to implant larger surface plate with lesser trauma and faster healing

Flexibility

Flexible tie-off capability in the silicone tube, which can be tight or loose to match complicated procedure

Unique Design

Unique fenestration design allow the growth of fibrous tissue to control bleb height and volumes, to secure plate in the place and minimize ocular motility disturbance

Radiographic Material

Barium impregnated silicone plate allows to identify with CT scan, MRI and plain skull film

Model No.	aadi 350
Surface Area	350 mm ²
Plate Length	31 mm
Tube Length	35 mm

Scan to view testimonials -
Study materials & Surgical videos



www.aurolab.com/glaucoma-shunt.asp

Information published in this catalogue is subject to change without notification.

I, Sivagangai Main Road, Veerapanjan, Madurai - 625 020, India
Tel : 91 452 309 6100 Fax : 91 452 244 6200
E.mail : info@aurolab.com www.aurolab.com

Issue 04 - 07/19



A R A V I N D E Y E C A R E S Y S T E M