

# Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)  
Absorbable, Synthetic, Monofilament

**Surgical Sutures**  
**Absorbable**



**Boz®**

## Definition

Short-term synthetic absorbable, closure of superficial wounds

## Material

Poly(glycolide (75%)- co-caprolactone (25%)) (PGCL)

## Coating

No coating

## Structure

Monofilament

## Colour

Undyed Beige (undyed) or Violet (dyed)

## Size

EP : 0.7 to 5

USP : 6/0 to 2

## Tensile Strength Retention

1. Week 45 % - 70 %
2. Week 20 %

## Mass Absorption

90 - 120 Days

## Sterilization Method



















Ethylene Oxide Gas

## Main Indication

Gynaecology  
Obstetrics  
Urology  
Plastic Surgery  
Odontology (oral mucosa)  
Skin closer - intracutaneous, subcutaneous.  
Reconstructive Surgery  
Gastrointestinal Surgery



















# Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)  
Absorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
8 mm							
1/4 Spatula 330 Micron			MK1398	6/0 0,7	45	12	Undyed
			MK1759	5/0 1,0	45	12	Undyed
10 mm							
3/8 Reverse Cutting 330 Micron			MK1494	6/0 0,7	45	12	Undyed
13 mm							
3/8 Round Body, Taper Point 330 Micron			MK1578	6/0 0,7	45	12	Undyed
3/8 Round Body, Taper Point 330 Micron			MK1905	5/0 1,0	50	12	Undyed
16 mm							
1/2 Round Body, Taper Point			MK1989	5/0 1,0	75	12	Undyed
			MK2503	4/0 1,5	75	12	Undyed
			MK3004	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MK2539	4/0 1,5	75	12	Undyed
3/8 Inside Cutting			MK2139	5/0 1,0	45	12	Undyed
			MK2557	4/0 1,5	75	12	Undyed
			MK3038	3/0 2,0	75	12	Undyed
18 mm							
3/8 Reverse Cutting			MK1704	6/0 0,7	45	12	Undyed
19 mm							
3/8 Reverse Cutting			MK2639	4/0 1,5	75	12	Undyed
			MK3108	3/0 2,0	75	12	Undyed













# Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)  
Absorbable, Synthetic, Monofilament

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body,Taper Point			MK2683	4/0 1,5	75	12	Undyed
			MK3141	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			MK2718	4/0 1,5	75	12	Undyed
			MK3167	3/0 2,0	75	12	Undyed
24 mm							
3/8 Inside Cutting			MK3236	3/0 2,0	75	12	Undyed
25 mm							
1/2 Round Body,Taper Point			MK2787	4/0 1,5	75	12	Undyed
			MK3251	3/0 2,0	75	12	Undyed
			MK3252	3/0 2,0	75	12	Violet
			MK3723	2/0 3,0	75	12	Undyed
3/8 Reverse Cutting			MK2805	4/0 1,5	75	12	Undyed
			MK3273	3/0 2,0	75	12	Undyed
3/8 Inside Cutting			MK3762	2/0 3,0	75	12	Undyed
26 mm							
1/2 Round Body,Taper Point			MK3312	3/0 2,0	75	12	Undyed
			MK3313	3/0 2,0	75	12	Violet
			MK3785	2/0 3,0	75	12	Undyed
1/2 Round Body,Taper Point			MK5846	2/0 3,0	75	12	Undyed
1/2 Reverse Cutting			MK3335	3/0 2,0	75	12	Undyed

# Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)  
Absorbable, Synthetic, Monofilament

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
26 mm						
<div>1/2</div> <div>Reverse Cutting</div> <div></div>	MK4229	0 3,5	75	12	Undyed	
	MK2851	4/0 1,5	75	12	Undyed	
<div>3/8</div> <div>Reverse Cutting</div> <div></div>	MK3345	3/0 2,0	75	12	Undyed	
	MK3821	2/0 3,0	75	12	Undyed	
30 mm						
<div>1/2</div> <div>Round Body,Taper Point</div> <div></div>	MK3368	3/0 2,0	75	12	Undyed	
	MK3369	3/0 2,0	75	12	Violet	
	MK3857	2/0 3,0	75	12	Undyed	
	MK3858	2/0 3,0	75	12	Violet	
	MK4243	0 3,5	75	12	Undyed	
	MK4607	1 4,0	75	12	Undyed	
<div>3/8</div> <div>Round Body,Taper Point</div> <div></div>	MK4263	0 3,5	75	12	Undyed	
<div>3/8</div> <div>Reverse Cutting</div> <div></div>	MK3890	2/0 3,0	75	12	Undyed	
	MK4272	0 3,5	75	12	Undyed	
	MK4628	1 4,0	75	12	Undyed	
35 mm						
<div>1/2</div> <div>Round Body,Taper Point</div> <div></div>	MK3915	2/0 3,0	75	12	Undyed	
	MK4299	0 3,5	75	12	Undyed	
	MK4643	1 4,0	75	12	Undyed	

# Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)  
Absorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>35 mm</b>					
<b>1/2</b> Reverse Cutting	<b>MK4658</b>	<b>1</b> 4,0	75	12	Undyed
<b>3/8</b> Reverse Cutting	<b>MK3932</b>	<b>2/0</b> 3,0	75	12	Undyed




<b>40 mm</b>					
<b>1/2</b> Round Body, Taper Point	<b>MK3981</b>	<b>2/0</b> 3,0	75	12	Undyed
	<b>MK4367</b>	<b>0</b> 3,5	75	12	Undyed
	<b>MK4374</b>	<b>0</b> 3,5	90	12	Undyed
	<b>MK4712</b>	<b>1</b> 4,0	75	12	Undyed
	<b>MK4713</b>	<b>1</b> 4,0	75	12	Violet
	<b>MK4718</b>	<b>1</b> 4,0	90	12	Undyed
	<b>MK5001</b>	<b>2</b> 5,0	75	12	Undyed
<b>1/2</b> Reverse Cutting	<b>MK4744</b>	<b>1</b> 4,0	75	12	Undyed
	<b>MK5018</b>	<b>2</b> 5,0	75	12	Undyed
<b>3/8</b> Reverse Cutting	<b>MK4407</b>	<b>0</b> 3,5	75	12	Undyed

<b>45 mm</b>					
<b>1/2</b> Round Body, Taper Point	<b>MK4771</b>	<b>1</b> 4,0	75	12	Undyed
	<b>MK5034</b>	<b>2</b> 5,0	75	12	Undyed

<b>50 mm</b>					
<b>1/2</b> Round Body, Taper Point	<b>MK4037</b>	<b>2/0</b> 3,0	75	12	Undyed
	<b>MK4811</b>	<b>1</b> 4,0	75	12	Undyed

# Monokaprol

Surgical Suture, Poly[glycolide(75%)-co-caprolactone(25%)](PGCL)  
Absorbable, Synthetic, Monofilament

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
<b>50 mm</b>					
<b>1/2</b> Round Body, Taper Point 	<b>MK4817</b>	<b>1</b> 4,0	90	12	Undyed
	<b>MK5065</b>	<b>2</b> 5,0	75	12	Undyed
<b>3/8</b> Reverse Cutting 	<b>MK4474</b>	<b>0</b> 3,5	75	12	Undyed
	<b>MK4840</b>	<b>1</b> 4,0	75	12	Undyed
<b>63 mm</b>					
<b>3/8</b> Blunt 	<b>MK4859</b>	<b>1</b> 4,0	100	12	Undyed
	<b>MK4857</b>	<b>1</b> 4,0	75	12	Undyed
Non-Needle Suture	Product No	USP EP	Length (cm)	Unit	Color
<b>150 cm</b>					
<b>MONOKAPROL</b> Non-Needle Suture	<b>MK2934</b>	<b>3/0</b> 2,0	150	12	Undyed
	<b>MK3510</b>	<b>2/0</b> 3,0	150	12	Undyed
	<b>MK4113</b>	<b>0</b> 3,5	150	12	Undyed
	<b>MK4531</b>	<b>1</b> 4,0	150	12	Undyed
	<b>MK4907</b>	<b>2</b> 5,0	150	12	Undyed