


DAHLHAUSEN®
Lo fut 86

Laryngeal Mask, single use

The sterile Laryngeal Mask, single-use, provides a respiration without intubation of the trachea. It seals the larynx to ensure safe respiration.

Material: Silicone, Polypropylene, Polycarbonat

Sterile: Yes

Art.no.	Application area	Max. filling volume (cuff)	Maximum body weight	Minimum body weight
4400060010	Newborn	4 ml	5 kg	<5 kg
4400060015	Small Child	7 ml	10 kg	5 kg
4400060020	Child	10 ml	20 kg	10 kg
4400060025	Child	14 ml	30 kg	21 kg
4400060030	Adolescent	20 ml	50 kg	31 kg
4400060040	Adult	30 ml	70 kg	51 kg
4400060050	Adult	40 ml	100 kg	71 kg

P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Str. 53
50996 Köln

