


TRIPLE BLOOD BAG



This system is used for separation of three components from whole blood. The triple system includes one primary bag with anticoagulant CPDA-1 solution BP and two empty satellite bags.

The triple SAGM system includes one primary bag with anticoagulant CPD Solution USP, SAGM solution in another bag and one empty satellite bag.

Sizes and Specification:



Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
TC-250W	250ml	2x200ml	16G	CPDA-1	35ml	2/50
TC-350W	350ml	2x300ml	16G	CPDA-1	49ml	2/50
TC-450W	450ml	2x300ml	16G	CPDA-1	63ml	2/50
TC-500W	500ml	2x300ml	16G	CPDA-1	70ml	2/50
TS-350W	350ml	2x300ml	16G	CPD+SAGM	49ml+80ml	2/50
TS-450W	450ml	2x300ml	16G	CPD+SAGM	63ml+ 102ml	2/50
TS-500W	500ml	2x300ml	16G	CPD+SAGM	70ml+ 114ml	2/50

Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
TC-250E	250ml	2x200ml	16G	CPDA-1	35ml	2/50
TC-350E	350ml	2x300ml	16G	CPDA-1	49ml	2/50
TC-450E	450ml	2x300ml	16G	CPDA-1	63ml	2/50
TC-500E	500ml	2x300ml	16G	CPDA-1	70ml	2/50
TS-350E	350ml	2x300ml	16G	CPD+SAGM	49ml+80ml	2/50
TS-450E	450ml	2x300ml	16G	CPD+SAGM	63ml+ 102ml	2/50
TS-500E	500ml	2x300ml	16G	CPD+SAGM	70ml+ 114ml	2/50

