

2- Suction Catheter (Thumb Control Connector)

Code No.: 2002



* Description:

⊙ Suction Catheter is used to provide access to respiratory tract to remove secretion with maximum comfort to the patient the connector regulates the vacuum by covering or partially covering a vent port with the thumb of the clinician who is doing the suctioning of the patient.

* Features:

- ⊙ For removal of Secretion from trachea and bronchial region.
- ⊙ Soft, frosted and kink resistant PVC tubing.
- ⊙ A traumatic, soft & rounded open tip with two lateral eyes.
- ⊙ Available with open end.
- ⊙ Color coded thumb control connector for easy identification of size.

*Options:

- ⊙ Tube Length: 52 cm.
- ⊙ Sizes: 5 FG to 20 FG.

* Material:

- ⊙ Medical Grade P.V.C. (DEHP free).
- ⊙ Master Batch.

<u>Code</u>	<u>Size</u>	<u>Color</u>
2002-5	5 FG	Grey
2002-6	6 FG	Light Green
2002-8	8 FG	Blue
2002-10	10 FG	Black
2002-12	12 FG	White
2002-14	14 FG	Green
2002-16	16 FG	Orange
2002-18	18 FG	Red
2002-20	20 FG	Yellow