

INJECTION NEEDLE


WHERE ERGONOMICS MEET RELIABILITY

MICRO-TECH's disposable injection needle enhances clinical practice due to its ergonomic handle and improved needle sharpness. The automatic "pen-style" handle allows single-handed use. For even more precise punctures, the needle has varying stiffness levels that can be gradually adjusted. When the needle is exposed, the quick-release


button snaps the needle back into the catheter to avoid damage to the working channel. The patented ratchet mechanism provides audible confirmation of a drawn needle, providing an effective exchange between the physician and nurse.

SPECIFIC CHARACTERISTICS

- Ergonomic handle with one-handed operation
- Patented ratchet mechanism with audible confirmation
- Adjustable stiffness levels
- Quick-release button for a safe needle retraction
- Enhanced needle sharpness for precise penetration
- Different needle diameters and lengths
- Easy attachment via a secure Luer-Lock connection









SPECIFICATIONS

| REF | Needle length mm | Needle Ø | Working channel Ø mm | Length mm |
|------------------------------------|------------------|---------------|----------------------|-----------|
| INJECTION NEEDLE - POLYPECTOMY/EMR | | | | |
| IN12-195232302 | 5 | 19 G = 1.1 mm | 2.8 | 2300 |
| IN12-225232302 | 5 | 22 G = 0.7 mm | 2.8 | 2300 |
| IN12-255232302 | 5 | 25 G = 0.5 mm | 2.8 | 2300 |
| INJECTION NEEDLE - ESD | | | | |
| IN12-224232302 | 4 | 22 G = 0.7 mm | 2.8 | 2300 |
| IN12-254232302 | 4 | 25 G = 0.5 mm | 2.8 | 2300 |

Packaging unit: 10 pieces