

BOND ELUT CERTIFY AND CERTIFY II

Part Numbers

Bond Elut Certify and Certify II are available in a wide variety of cartridge volumes to accommodate both manual and automated drug extraction methods. Sorbent masses from 50 mg to 1 g are available for both standard drug testing methods and pharmaceutical research.

Bond Elut Certify and Certify II are also available in 96-well plate formats. These formats include our 1 mL 96-well plate, 2 mL 96-square well plate, and the VersaPlate format. The VersaPlate format incorporates 96 individual SPE cartridges that can be assembled into a 96-well array, using a reusable holder, or the VersaPlate tubes function as standalone SPE cartridges for vacuum or positive pressure processing. The VersaPlate array format minimizes waste, ensures high throughput, and is ideal for scalability in your methods.

Bond Elut Certify

Description	Unit	40 μ m Particle Size	120 μ m Particle Size
Large Reservoir Capacity (LRC) Cartridges			
130 mg, 10 mL	50/pk	12113050	14113050
130 mg, 10 mL	500/pk	52113050	14113055
200 mg, 10 mL	50/pk	12113054	14113054
200 mg, 10 mL	500/pk	52113051	
300 mg, 10 mL	50/pk	12113052	14113052
Straight Barrel Cartridges			
50 mg, 3 mL	50/pk	12105030	
130 mg, 1 mL	100/pk	12102083	14102083
130 mg, 3 mL	50/pk	12102051	14102051
130 mg, 3 mL	500/pk	52102051	
130 mg, 3 mL ASPEC	500/pk	22102084	
130 mg, 3 mL tabless	50/pk	12102051T	
130 mg, 6 mL	30/pk	12256146	
130 mg, 6 mL tabless	500/pk	12256146TJ	
200 mg, 3 mL	50/pk	12102145	
200 mg, 6 mL	30/pk	12256145	
300 mg, 3 mL	50/pk	12102081	
300 mg, 3 mL	500/pk	52102081	
300 mg, 3 mL ASPEC	500/pk	22102081	
300 mg, 3 mL tabless	50/pk	12102081T	14102081T
300 mg, 6 mL	30/pk	12102082	
500 mg, 6 mL	30/pk	12102093	14102093
1 g, 6 mL	30/pk	12102085	14102085

TIP: To meet your lab's high-throughput needs, choose cost-effective bulk packaging for many of our products. Visit [agilent.com/chem/bigpacks](https://www.agilent.com/chem/bigpacks)