

HOME ABOUT US

CATALOG OEM-OBL-PL PRODUCTS V

**NEWS FROM US** 

# LUBRICANT GEL WITH LIDOCAIN

HOME > MEDICAL > LUBRICANT GEL WITH LIDOCAIN



### ONEGEL

# LUBRICANT GEL WITH LIDOCAIN

# **USAGE AREAS**

ONEGEL is a sterile, water-soluble lubricating gel used to lubricate catheters and other medical devices during urethral, colon and rectal applications such as urethral catheterization, endoscopy and cystoscopy.

It can be also used for rectal and colonic applications as a lubricant gel.

The most important function of the gel is to form a lubricating layer between the catheter or other medical device and the urethral, rectal and colonic mucosa.

# PRODUCT PROPERTIES

Disposable.

Soluable in water.

Easy to use thanks to the syringe and accordion tube.

Eliminates the risk of damage to tissue.

Contains antiseptic Chlorhexidine Gluconate to eliminate the risk of infection.

Lidocaine Hydrochloride is a local anesthetic agent used to prevent pain.

It provides a painless catheterization with its anesthetic effect.

WhatsApp



Thanks to the lubricating feature provided by hydroxyethylcellulose, the patient relaxes and minimizes the spasticity and



HOME ABOUT US CATALOG OEM-OBL-PL PRODUCTS V NEWS FROM US C

Remove the stopper at the end of the syringe/accordion tube.

Pour a drop of gel to facilitate application.

After placing the syringe tip on the area to be lubricated, gently squeeze some ONEGEL Lidocaine Lubricating Gel by pressing the plunger of the syringe/accordion tube.

## **COMPOSITION**

Pure Water

Propylene Glycol, Hydroxyethylcellulose (Lubricant)

Lidocaine Hydrochloride (Local Anesthetic)

Chlorhexidine Gluconate - 20% Concentration (Antiseptic)

Methyl Hydroxybenzoate (Preservative)

Propyl Hydroxybenzoate (Preservative)

#### **SYMBOLS**

Contains 2% Lidocaine.

STERILE.

Antiseptic.

Anesthetic

It is a Class III Medical Device within the scope of 93 / 42 / EEC Medical Device Regulation.

#### TABLE:


Product Name ONEGEL 6ML ONEGEL 11ML ONEGEL 12,5GR

Product Reference Code KAF G27-6 KAF G27-11 KAF G27-12,5

Packaging Type Syringe Syringe Accordion Tube

GMDN 37717 37717 37717

Packaging 25 x 6ml / Box 25 x 11ml / Box 25 x 12,5gr / Box





HOME ABOUT US CATALOG OEM-OBL-PL PRODUCTS  $\lor$  NEWS FROM US (

WhatsApp QGetButton