

Huaian Helen Medical Instrument Co., Ltd

Address: Group 3, Zhengtai Village Matou Town, Huaiyin District Huaian, 223345 Jiangsu, China, Tel/Fax: 0517-84266088

Ce Approved Medical Supply Single-Use Devices Safety Pressure Activated Blood Lancet



Disposable blood collection needles are divided into four types: twist blood lancet, flat blood lancet, safety pressure blood lancet and stainless steel lancet. The twist lancet and flat lancet are composed of a steel needle and a needle. The safety pressure blood lancet lancet is composed of a needle body, a spring and a plastic part. The needle body is composed of a steel needle and a needle. The twist blood lancet, flat blood lancet, safety pressure blood lancet and stainless steel lancet follow the diameter of the steel needle can be divided into 21G, 23G, 26G, 28G, 30G. The slice lancet is made of stainless steel, and is divided into two specifications: large and small according to different sizes. The product is sterilized by gamma radiation, and the product should be sterile.

- 1. Unique single-use design maximizes protection and the needle is completely shielded before and after use, thus minimizing the risk of needle injury and cross-contamination between different people.
- 2. Sterility: Gamma ray, safety cover to protect the needle, built-in needle protection system.
- 3. Minimize tissue damage: Linear needle tracking can eliminate needle vibration, and the release button remains pressed after use to ensure one-time operation.
- 4. The high quality tri-bevel penetration method minimizes patient pain and clinician error.

Model	Sterile	Material	Depth	Gauge	Packing	Ctn size (cm)	GW (kg)
XY21G XY26G XY28G XY30G	Gamma	Stainless Steel	1.8- 2.4mm	21g, 26g, 28g, 30g	100pcs/box, 20box/ctn	33*28.4*15.5	6.6
XY23G				23g			