

elastopor STERIL non-woven dressing with absorbent pad, self-adhesive, sterile

TD-09

att.no.3 pt.2 subpt.2.3

edition II

date: 2020-01-10

change of day: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- Treating wounds with a small and medium amount of exudate
- Especially recommended for covering surgical wounds characterized by straight cuts and the application of fewer or more skin sutures in a straight line

- Self-adhesive dressing with absorbent pad
- Made of hydrophobic non-woven material
- Micropores in the non-woven material structure provide an adequate rate of moisture vapour permeability so correct gas exchange between the dressing and skin is guaranteed
- Centrally placed absorbent pad absorbs medium wound exudate
- Absorbent pad coated with polyethylene mesh protects against external factors' influence and does not stick to the wound as well as facilitates removal of the dressing
- Hypoallergenic
- Elastic comfortable for the patient, adjusts to the body shape and guarantees freedom of movement
- Easy to attach to the skin
- Rounded edges prevent accidental removal of the dressing
- Covered with hypoallergenic acrylic glue without the content of zinc oxide, rubber, latex, evenly applied over the entire surface, with very good moisture resistance ensuring good stabilization of the dressing at the application site
- Silicone paper backing facilitates precise, painless and efficient application
- EO sterilized



SIZE AND REFERENCE						
CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	INTERMEDIATE PACKAGING BOX	BULK PACKAGING		
801001	5 cm x 7,2 cm	1 pc.	100 pcs	24 x 100 pcs		
801012	6 cm x 10 cm	1 pc.	50 pcs	45 x 50 pcs		
801002	8 cm x 10 cm	1 pc.	30 pcs	36 x 30 pcs		
801013	10 cm x 10 cm	1 pc.	30 pcs	24 x 30 pcs		
801021	10 cm x 12 cm	1 pc.	30 pcs	30 x 30 pcs		
801014	10 cm x 15 cm	1 pc.	30 pcs	24 x 30 pcs		
801003	6 cm x 15 cm	1 pc.	50 pcs	24 x 50 pcs		
801004	8 cm x 15 cm	1 pc.	30 pcs	24 x 30 pcs		
801005	8 cm x 20 cm	1 pc.	30 pcs	18 x 30 pcs		
801015	10 cm x 20 cm	1 pc.	25 pcs	18 x 25 pcs		
801016	10 cm x 25 cm	1 pc.	25 pcs	18 x 25 pcs		
801017	10 cm x 30 cm	1 pc.	25 pcs	12 x 25 pcs		
801018	10 cm x 35 cm	1 pc.	25 pcs	12 x 25 pcs		



elastopor STERIL non-woven dressing with absorbent pad, self-adhesive, sterile

TD-09

att.no.3 pt.2 subpt.2.3

edition I

date: 2021-10-29

change of day: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- treating wounds with a small and medium amount of exudate
- especially recommended for covering surgical wounds characterized by straight cuts and the application of fewer or more skin sutures in a straight line

- self-adhesive dressing with absorbent pad
- made of hydrophobic non-woven material
- micropores in the non-woven material structure provide an adequate rate of moisture vapour permeability so correct gas exchange between the dressing and skin is guaranteed
- centrally placed absorbent pad absorbs medium wound exudate
- absorbent pad coated with polyethylene mesh protects against external factors' influence and does not stick to the wound as well as facilitates removal of the dressing
- hypoallergenic
- elastic comfortable for the patient, adjusts to the body shape and quarantees freedom of movement
- easy to attach to the skin
- rounded edges prevent accidental removal of the dressing
- covered with hypoallergenic acrylic glue without the content of zinc oxide, rubber, latex, evenly applied over the entire surface, with very good moisture resistance ensuring good stabilization of the dressing at the application site
- silicone paper backing facilitates precise, painless and efficient application
- EO sterilized



SIZE AND REFERENCE					
CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	INTERMEDIATE PACKAGING BOX	BULK PACKAGING CARTON	
801001E	5 cm x 7,2 cm	1 pc.	100 pcs	24 x 100 pcs	
801012E	6 cm x 10 cm	1 pc.	50 pcs	45 x 50 pcs	
801002E	8 cm x 10 cm	1 pc.	30 pcs	36 x 30 pcs	
801013E	10 cm x 10 cm	1 pc.	30 pcs	24 x 30 pcs	
801021E	10 cm x 12 cm	1 pc.	30 pcs	30 x 30 pcs	
801014E	10 cm x 15 cm	1 pc.	30 pcs	24 x 30 pcs	
801003E	6 cm x 15 cm	1 pc.	50 pcs	24 x 50 pcs	
801004E	8 cm x 15 cm	1 pc.	30 pcs	24 x 30 pcs	
801005E	8 cm x 20 cm	1 pc.	30 pcs	18 x 30 pcs	
801015E	10 cm x 20 cm	1 pc.	25 pcs	18 x 25 pcs	
801016E	10 cm x 25 cm	1 pc.	25 pcs	18 x 25 pcs	
801017E	10 cm x 30 cm	1 pc.	25 pcs	12 x 25 pcs	
801018E	10 cm x 35 cm	1 pc.	25 pcs	12 x 25 pcs	



elastopor EYE eye dressing, non-woven, with absorbent pad, self-adhesive, sterile

TD-16

edition I

att.no.3 pt.2 subpt.2.3

date: 2016-06-03

change of day: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- Covering of an eye and its area after ophthalmic procedures and surgeries
- Sight correction in strabismus therapy

- Self-adhesive eye dressing with absorbent made of hydrophobic non-woven material
- Covered with acrylic glue
- Micropores in the non-woven material structure provide an adequate rate of vapour permeability so correct gas exchange between the dressing and skin is guaranteed
- Centrally situated absorbent pad absorbs moderate amounts of exudates and protects against external influences
- Absorbent pad with polyethylene mesh does not stick to the
- Elastic comfortable for the patient, adjusts to the body shape and guarantees freedom of movement
- Rounded edges prevent accidental removal of the dressing
- Silicone paper backing facilitates precise, painless and efficient application
- Available in various sizes and colours of non-woven and absorbent pad
- Hypoallergenic
- Sterile



SIZE AND REFERENCE						
CATALOGUE NUMBER	DIMENSIONS	SHAPE	COLOUR	INDIVIDUAL PACKAGING PAPER-PAPER	INTERMEDIATE PACKAGING BOX	BULK PACKAGING
803010	5,8 cm x 8,3cm		absorbent pad – white	1 pc.	50 pcs	24 x 50 pcs
803011	6,5 cm x 9,5 cm		non-woven material - white	1 pc.	50 pcs	24 x 50 pcs
803212	5 cm x 7,5 cm		absorbent pad - white non-woven material - skin colour	1 pc.	50 pcs	24 x 50 pcs



elastoLUMENAL S Eye pad, superabsorbent, multi-layer, sterile

att.no.3 pt.2 subpt.2.3

edition I

date: 2016-06-03

change of day: 2021-12-22

TD-16

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- Covering an eye and its area after ophthalmic procedures and surgeries
- Sight correction in strabismus therapy

- Composed of external non-woven layers and internal multi-layer cellulose pulp
- Needs to be fixed with adhesive tape or in different way, e.g. with tubular net dressing
- Hypoallergenic
- Sterile



		SIZE AND REFERENCE		
CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	INTERMEDIATE PACKAGING BOX	BULK PACKAGING
803310	5,5 cm x 7,5 cm	1 pc.	50 x 1 pcs	24 x 50 pcs



elastopor IV IV cannula dressing, non-woven, self-adhesive, sterile

TD-09

att.no.3 pt.2 subpt.2.3

edition II

date: 2020-01-10

change of day: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

 Stabilization, protection and fixation of intravenous cannulas with a vertical port

- Self-adhesive dressing having a notch for the vertical port of the intravenous cannula
- Made of hydrophobic non-woven material
- Micropores in the non-woven material structure provide an adequate rate of moisture vapour permeability so correct gas exchange between the dressing and skin is guaranteed
- Centrally situated absorbent pad protects the nonwoven material from sticking to the injection site
- Non-woven lining should be placed under cannula's wings, which prevents it from causing damage to the skin
- Hypoallergenic
- Elastic adjusts to the body shape and guarantees freedom of movement
- Easy to attach to the skin
- Rounded edges prevent accidental removal of the dressing
- Coated with a hypoallergenic acrylic adhesive without zinc oxide, applied to the nonwoven fabric by a continuous method
- Silicone paper backing facilitates precise, painless and efficient application
- Graphical instructions for use on individual and intermediate packaging
- EO sterilized
- Shelf life: 5 years



		SIZE AND REFERENCE		
CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	INTERMEDIATE PACKAGING BOX	BULK PACKAGING
811013	6 cm x 8 cm	1 pc.	100 pcs	20 x 100 pcs
811012	5,1 cm x 7,6 cm	1 pc.	100 pcs	20 x 100 pcs
811023	6 cm x 8 cm	1 pc.	50 pcs	40 x 50 pcs
811022	5,1 cm x 7,6 cm	1 pc.	50 pcs	40 x 50 pcs



Adhesive silk tape SILKplast

TD-06-I.2.c-5.2

edition I

date: 2021-05-26

amendment of: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS:

ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com

CLASSIFICATION

class I medical device compliant with Medical Devices Regulation (EU) 2017/745 (MDR)

INTENDED PURPOSE

 device intended for stabilisation, protection and fixation of drains, catheters, medical tubes not requiring specialised fixation at their place of application as well as primary dressings (contact), non-adhesive, applied directly to wounds and skin damages

- made of artificial silk
- white colour
- wound on a roll
- covered with an acrylic adhesive applied by the continuant method
- can be divided lengthwise and crosswise without the use of scissors
- serrated edge finishing
- good air and water vapour permeability
- good adhesion
- thin, flexible, conforms to the body contours, guarantees freedom of movement
- hypoallergenic
- single-use
- non-sterile
- shelf life: 3 years



CATALOGUE NUMBER	DIMENSIONS	INTERMEDIATE PACKAGING (SALES) BOX	BULK PACKAGING CARTON
SS12505W	1,25 cm x 5 m	24 x 1 pc.	30 x 24 pcs
SS12509W	1,25 cm x 9,14 m	24 x 1 pc.	30 x 24 pcs
SS25005W	2,5 cm x 5 m	12 x 1 pc.	30 x 12 pcs
SS25009W	2,5 cm x 9,14 m	12 x 1 pc.	30 x 12 pcs
SS50005W	5 cm x 5 m	6 x 1 pc.	30 x 6 pcs
SS50009W	5 cm x 9,14 m	6 x 1 pc.	30 x 6 pcs