


### Teveik Part Number Code Rule:



a.—— Product Line Code

1- Indicates Reusable SpO2 Sensor    2- Indicates Disposable SpO2 Sensor

b.—— Brand Code

For Example: 153-Indicates Nellcor    164-Indicates Philips

c.—— Difference Code For Different Machine Model From Same Brand

With Letter: A; B; C....

d.—— Sensor End Code ( Refer to next page)

1: Adult finger clip	4: Child soft tip	7: multi-site Y type	A: multi-site clip
2: Child finger clip	5: Infant soft tip	8: adult ear clip	B: adult Forehead
3: adult soft tip	6: Neonate silicone wrap	9: veterinary tongue clip	C: Y type with foam

e.—— Cable length


3: L=3m

f.—— Difference code for cable


A-Indicates blue TPU single cable	B-Indicates blue TPU double cable
C-Indicates grey TPU single cable	D-Indicates grey TPU double cable

For example:

1147E13A : Mindray adult finger clip SpO2 sensor, redel 7pin,10ft



### Teveik Part Number Code Rule:



a. — Product Line Code

1- Indicates Reusable SpO2 Sensor 2- Indicates Disposable SpO2 Sensor

b. — Brand Code

153-Indicates Nellcor 164-Indicates Philips

c. — Difference Code For Different Machine Model From Same Brand

With Letter: A; B; C....

d. — Patient Sizes Code

1: Adult	3: Infant
2: Child	4: Neonate

e. — Material Code For Sensor End









1: Medaplast	3: Microfoam
2: Transpore	4: velcro

f. — Difference Code








With Letter: A; B; C....

For example:

2153E12A: Nellcor adult disposable SpO2 sensor,0.45m,DB7pin,transpore

Sensor Type	Picture	Patient Sizes	Application Site
Medaplast		Adult: > 40kg Pediatric: 10-40kg Infant: 3-15kg Neonate: 1-4kg,	
Transpore		Adult: > 40kg Pediatric: 10-40kg Infant: 3-15kg Neonate: 1-4kg,	
Microfoam		Adult: > 40kg Pediatric: 10-40kg Infant: 3-15kg Neonate: 1-4kg,	
velcro		Adult: > 40kg Neonate: 1-4kg,	



Compatible Brand	Compatible Model	Connector Type	OEM P/N	Teveik P/N	Spec
BCI	BCI: 3100 3101 3300 3301 3302 3303 3304 3401 6100 9100, MiniTorr, MiniTorr Plus, Drager:Oxipac 2000/2500, Spacelabs Burdick:PRO2,100 Aemsrong:AD- 900,1000,2000,3000;Corpluse:08/16; Digicare:Digimax5500; Invacare;all,Lynn medical:all;Medlab:POX10;ODAM:AL L	 DB 7pin		2027C1*A	adult
				2027C2*A	Pediatric
				2027C3*A	Infant
				2027C4*A	neonate
CSI	503DX,503Spot,504DX,506,506DX	 DB 6pin		2060E1*A	adult
				2060E2*A	Pediatric
				2060E3*A	Infant
				2060E4*A	neonate
CSI		 DB 6pin		2060D1*A	adult
				2060D2*A	Pediatric
				2060D3*A	Infant
				2060D4*A	neonate
Datascope		 DB 9pin		2065D1*A	adult
				2065D2*A	Pediatric
				2065D3*A	Infant
				2065D4*A	neonate
Datex	AS/3, CS/3,S/5,CH/S, CH/2S, CH- RS, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Capnomac Ultima, Light Monitor, M-Estpr, Oscar II, OSP-200, Satlite Trans, Satlite Trans Plus, ULTSVI, ULTSVO, ULTSIO, ULTS00.	 DB 9pin		2066D1*A	adult
				2066D2*A	Pediatric
				2066D3*A	Infant
				2066D4*A	neonate
Dolphin		 DB 9pin		2073C1*A	adult
				2073C2*A	Pediatric
				2073C3*A	Infant
				2073C4*A	neonate
Ohmeda	Ohmeda:Tuffsat,3775,37xx series,3800,3900,3900P,4700,5250R GM,M-OSAT module, N-XOSAT option for Cardiocap@/5,,RASCALII, Mennen: Horizon XL800 (Ohmeda spo2);GE : 7000 SERIES,COROMETRICS 511,TRAM A.R.,TRAM MODULAR X00SL SERIES,as 200SL,400SL,600SL,800A,800SL	 DB 9pin		2067D1*A	adult
				2067D2*A	Pediatric
				2067D3*A	Infant
				2067D4*A	neonate