

Sales only to commercial customers. No sales to consumers within the meaning of §13 BGB.



Zoom in

One-bottle Complete Set, with Heimlich-valve

The Thorax-drainage System is used for collection of air and/or liquid from the pleural cavity. Set Composition: One bottle 1.300 ml, patient tube 150 cm –patient-connector with puncture membrane, vacuum tube 100 cm with stepped-funnel, Heimlich valve.

Set content:

One bottle; one Heimlich-valve; one 150 cm patient tube; one 100 cm vacuum tube; connector with puncture membran

Sterile:

Yes

Country availability:

Please note: not all products are available in all countries. Feel free to contact our

[read more](#)

Article no.

0410111150

minimal quantity

15