

# MINDRAY TV80 NEONATAL TRANSPORT VENTILATOR

## Accessory and Supply List - Without Commercial Information

<b>Equipment</b>	Mindray TV80 Transport Ventilator
<b>Configuration</b>	Neonatal transport ventilator configuration
<b>Purpose of this document</b>	List of accessories and supply items included with the proposed neonatal ventilator configuration, excluding prices and commercial conditions.

### 1. Main equipment and included configuration

Item code	Description	Quantity	Included configuration / remarks
ET20B-PA00009	TV80 Ventilator (VS + APRV + AMV + Neo)	1 unit	TV80 main unit, Spanish language configuration, EU power cord, AC adapter, O2 hose 3 m Fra/NIST, NIST connector, non-consuming oxygen sensor, module slot, Bluetooth, TV ventilator system software, neonatal function, neonatal disposable circuit kits and neonatal demo lung.

### 2. Neonatal disposable patient connection sets

Item code	Description	Quantity	Included configuration / remarks
115-096082-00	Disposable neonatal circuit, dual limb, with flow sensor for TV80	20 sets	Each set includes one disposable neonatal dual-limb circuit and one disposable neonatal flow sensor, differential pressure type. Intended as neonatal patient connection set for Mindray TV80.

### 3. Summary statement for tender compliance

The proposed Mindray TV80 neonatal transport ventilator configuration includes 20 neonatal patient connection sets. Each set is composed of a disposable neonatal dual-limb circuit and a disposable neonatal flow sensor for TV80.