



## TECHNICAL DATA SHEET

Document No	DD.8.2.6.2.02
Validity	01.06.2021
Rev No / Date	-
Page No	1 / 1

**ITEM** : 3 PARTS 5 ML SYRINGE with 22G x 1 1/4" NEEDLE

**DESCRIPTION** : Sterile hypodermic syringe for single use that is made of plastic material, pyrogen free, non-toxic, sterilized with ethylene oxide and used for aspiration of liquids or used for injection of liquid after filled

### PROPERTIES

#### PHYSICAL PROPERTIES

	STANDARD ( TEST METHOD)	UNITE	RESULT
Material	TS EN ISO 7886-1	-	Medical
Capacity	TS EN ISO 7886-1	ml	5 (± % 5)
Death Space	TS EN ISO 7886-1	ml	Max.0,075
Total graduated capacity	TS EN ISO 7886-1	ml	5 ml
Dimensions	TS EN ISO 7886-1	mm	Ø inside = 12,40 Ø outside = 13,10 l barrel = 62,00 l plunger = 73,50 l plunger + gasket = 75,70
Air Leakage during aspiration	TS EN ISO 7886-1	kPa/sn	No leakage at 88 kPa vacuum at 60 second
Liquid leakage under compression	TS EN ISO 7886-1	kPa/sn	No leakage at 300 kPa pressure at 30 second
Conic fittings (centric)	TS 3521-1 EN 20594-1	%	6% (Luer)
Needle size	TS 4002 EN ISO 7864	mm	0.70 x 30
Needle penetration force	TS 4002 EN ISO 7864	N	0.48-0.53
Bond between needle and tube	TS 4002 EN ISO 7864	N	Does not broken by 40N
Cylinder-piston difference	TS EN ISO 7886-1 (Article 12)	mm	15 mm
Scale length	TS EN ISO 7886-1 (Article 10.3)	mm	5 ml : 41,4 mm
Scale interval	TS EN ISO 7886-1 (Article 10.2)	ml	0,2 ml
Increment between graduation lines to be numbered	TS EN ISO 7886-1 (Article 10.2)	ml	1 ml
Lubricant	TS EN ISO 7886-1	gr/cm <sup>2</sup>	0,25 polydimethyl siloxane
Weight	AldarMed	gr	4,50-4,65

#### CHEMICAL PROPERTIES

Neutrality	TS EN ISO 7886-1	pH	Difference less than 1
Pb, Zn, Sn, Fe	TS EN ISO 7886-1	mg/lt	Totally less than 5
Cd	TS EN ISO 7886-1	mg/lt	Less than 0,1
Ethylene oxide residue	TS 9897	mg/lt	Less than 5

#### BIOLOGICAL PROPERTIES

Cytotoxicity	EN ISO 10993		Not cytotoxic
Toxicity	EN ISO 10993		Does not cause acute systemic toxicity
Pyrogenity	EN ISO 10993-11. European Pharmacopoeia 7 <sup>th</sup>		Apyrogen
Sterility	EN ISO 556, TS 10409/October 1992		Sterile
Hemolytic effect	EN ISO 10993-11, TS 6459/January		Not hemolytical
Sensitization	EN ISO 10993-10/March 1998		No sensitivity
Irritation	EN ISO 10993-10/March 1998		No irritation effect

#### PACKAGING PROPERTIES AND HANDLING

Primary Packing	TS EN ISO 7886-1	-	Transparent film + Blister Packing from medical grade paper and statements expressed in the standards
Secondary Packing	TS EN ISO 7886-1	-	inner and out surface testliner, fluting wave, offset printed box and statements expressed in the standards
Transport-Storage-Packing	AldarMed	-	6 inner boxes are put in outer boxes or wrapped container
Storage Condition	AldarMed	-	Kept out of direct sun light and stored in clean and dry storage areas at max. 45 °C

	<b>PREPARED BY</b>	<b>APPROVED BY</b>
<b>NAME &amp; SURNAME</b>	Sümeyye Ersoy	Seydi Nur Talay
<b>DUTY</b>	Quality Control Responsible	Quality Manager

