LeLun, Life Guard!





工厂地址:江苏省常州市横山桥镇

Factory Add: Hengshanqiao Town, Changzhou, Jiangsu, China办公地址:中国常州关河东路66号九洲环宇商务广场1808室

Office Add: Room 1808, 66 East Guanhe Road, ECO Business Plaza Changzhou, China

Tel: +86-519-88168398 / 88602288 Fax: +86-519-88168398 / 88603638

E-mail: info@lelun.com Website: www.lelun.com



LELUN MEDICAL APPLIANCES 乐论医疗

CHANGZHOU MEDICAL APPLIANCES GENERAL FACTORY CO.,LTD.

About US

Changzhou Medical Appliances General Factory Co., Ltd was built in 1988, it is a modern factory specialized in producing disposable medical appliances in China. The traffic is convenient and close to Shanghai.

The area of the factory is 50000m², the area of purifying workshop is 15000m², and fixed assets are about 5,000,000USD. Our main products are Disposable syringe, Safety syringe, Auto-disable syringe, Disposable infusion set, Hernia mesh, Medical stapler, Disposable blood transfusion sets, Urine bag, IV cannula, Oxygen mask, Examination glove, Surgical glove, Vaginal Speculum etc. Now we can manufacture more than 100,000,000 sets per year.

Through 30 years hard work, our hardware and software have up to the first-class in the same factory in China. And we won the title of "High quality product" by Jiangsu Province in 1991. In 2000, our factory passed through the ISO13485 and got the certificate of CE. In 2005, We had the license of exporting by ourselves. In 2011, we passed the GMP by Chinese SFDA.

Now our products are not only sold to Chinese market, but also sold to more than 60 countries like Germany, USA, Canada, UK, Italy, Austria, Romania, France, Bulgaria, Ghana, Morocco, Turkey, Colombia, Mexico, Bangladesh, Peru, Chile, Sri Lanka, Azerbaijan, South Africa, Libya, Sudan, Yemen, Pakistan, Togo and so on.

Our factory persists in the principle of "first quality, high credit, and superior service". Sincerely welcome customers all over the world to talk business and cooperate with us.





















Production Workshop and Equipment

ADVANCED EQUIPMENTS IS CONSOLIDATE THE BASIS OF QUALITY











CONTENTS

Disposable Syringe Series page 01~04 page 05~06 Oral / Enteral Syringe Series page 07~10 **Infusion Set Series Needle Series** page 13~14 **Urine / Stool Container Series Mask Series** page 15 **Medical Glove Series** page 16Hernia Mesh Series page 16 page 17Components



Disposable Three-Part Syringe Luer Slip Without Needle

• Item Code: CMSSxx



Disposable Three-Part Syringe Luer Slip With Needle

- Item Code: CMSSxxGx
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml



Disposable Three-Part Syringe Luer Lock Without Needle

- Item Code: CMSLxx
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml, 100ml, 200ml

Disposable Three-Part Syringe Luer Lock With Needle

- Item Code: CMSLxxGx
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml, 100ml, 200ml



Two-Part Syringe Without Needle

- Item Code: CMSTxx
- Size: 2ml, 5ml, 10ml, 20ml

Two-Part Syringe With Needle

- Item Code: CMSTxxGx
- Size: 2ml, 5ml, 10ml, 20ml



Insulin Syringe

- Item Code: CMSI
- Size: 0.3ml, 0.5ml, 1ml(U-100/U-40)
- Fixed needle size: 27g, 28g, 29g, 30g, 31g, 32g
- Fixed needle length: 6mm, 8mm, 13mm

Irrigation Syringe

- Item Code: CMSC
- Size: 20ml, 30ml, 50ml, 60ml, 100ml, 120ml, 150ml, 200ml, 250ml, 300ml, 350ml, 500ml
- 50ml/60ml with bulb
- 100ml with ring



Pump Syringe

- Item Code: CMSP
- Size: 50/60ml





Safety Syringe with retractable needle (Manually)

- Item Code: CMSA
- Size: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml
- Simple hypodermic needle available
- Needle be retractable into the barrel and locked after use





Auto Disable Syringe


- Item Code: CMSAD
- Size: 0.05ml, 0.1ml, 0.3ml, 0.5ml, 1ml, 2/3ml, 5ml, 10ml
- Simple hypodermic needle available
- Plunger be locked in the barrel after use



Veterinary Syringe

- Item Code: CMSV
- Size: 5ml, 8ml, 10ml, 13ml, 20ml, 30ml, 60ml
- Main Usage: 1. Cow mastitis syringe
 - 2. Goat & Sheep Wormer syringe

 - 3. Oral paste syringe for horse 4. Feeding syringe for calves to supplement
 - 5. Insecticide gel syringe





ORAL / ENTERAL SYRINGE SERIES

Oral / Enteral Syringe with Tip Cap

- Item Code: CMOSC
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- Non-luer tip, can not fit with hypodermic needles
- Latex or latex free, medical grade PP
- Transparent barrel, colored plunger
- With tip cap for racap after use



Oral / Enteral Syringe with Adaptor

- Item Code: CMOSA
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- Non-luer tip, can not fit with hypodermic needles
- Latex or latex free, medical grade PP
- Transparent barrel, colored plunger
- With transparent adaptor for bottle, fitting the bottle neck ID 14mm~24mm



Oral/Enteral Syringe Fernale Reverse Luer Lock

- Item Code: CMOSLL
- Size: 1ml, 3ml, 5ml, 10ml, 20ml, 30ml, 60ml
- Fearture: Enfit syringe, High transparency, accurate and clear scale marking



Oral Syringe (Luer Slip Tip) with Dispensing Tube

- Item Code: CMOSDT(LS)
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- Luer tip
- With protective cap and dispensing tube



Oral Syringe (Two Parts) with Adaptor

- Item Code: CMOSA-1
- Size: 2/3ml, 5ml, 10ml, 20ml
- Luer tip
- With adaptor fitting the bottle neck ID 18mm~22mm





- Type: Luer slip or luer lock
- Material: Medical grade PVC
- Specification: 1. Drip chamber with air vent
 - 2. Medical grade PVC tube, 120~180mm
 - 3. Rubber drug-adding port
 - 4. Polybag or blister packing
 - 5. With hypodermic needle or with scalp vein set
 - 6. Drip chamber: 48mm or 54mm, 20 drips / 1ml for adult; 60 drips / 1ml for pediatric
 - 7. Gravity use only



Disposable Infusion Set with Y-connector


- Item Code: CMIF-2
- Type: Luer slip or luer lock
- Material: Medical grade PVC
- Specification: 1. Drip chamber with air vent
 - 2. Medical grade PVC tube, 120~180mm
 - 3. Y drug-adding port
 - 4. Polybag or blister packing
 - 5. With hypodermic needle or with scalp vein set
 - 6. Drip chamber: 48mm or 54mm, 20 drips / 1ml for adult; 60 drips / 1ml for pediatric
 - 7. Gravity use only



Burette Infusion Set

- Item Code: CMIFB
- Size: 100ml, 150ml







Infusion Set for Pump

- Item Code: CMIP
- Material: medical grade PVC
- Specification: 1. Small / big chamber with air vent
 - 2. Medical grade PVC tube, 120~180mm
 - 3. Polybag or blister packing
 - 4. Drip chamber: 48mm or 54mm, 20 drips / 1ml

With 200 micron filter and 18 G hypodermic needle

5. It can be used with the pump



Size		Color		Length					
OD(mm)	Gauge	Color & Code		1/2" (13mm)	5/8" (15mm)	3/4" (20mm)	1" (25mm)	1 1/4" (32mm)	1 1/2" (38mm)
0.3	30G		Light yellow	√					
0.33	29G		Red	√					
0.36	28G		Blue-green	√					
0.4	27G		Grey	$\sqrt{}$					
0.45	26G		Brown	√	√				
0.5	25G		Orange		√	√	√		
0.55	24G		Purple			√	√		
0.6	23G		Blue				√	√	√
0.7	22G		Black				√	√	√
0.8	21G		Green				√	√	√
0.9	20G		Yellow				√		
1.1	19G		Beige				√		
1.2	18G		Pink				√		√
1.6	16G		White						√



Scalp Vein Set

- Item code: CMIV
- Size: 16G/17G/18G/19G/20G/21G/22G/ 23G/24G/25G/26G/27G
- Luer slip or luer lock
- Butterfly wings can make better handling and attachment to skin



I.V. Cannula

- Item Code: CMCI-2
- Size: 14G/16G/18G/20G/22G/24G/26G

Spinal Needle

- Item Code: CMSN
- Quincke Bevel & Pencil Point
- Size: 18G/19G/20G/21G/22G/23G/24G/25G/26G/27G
- Specification: 90mm





Item Code: CMBC

• Size: 20G/21G/22G/23G

Item Code: CMSNS

11/www.lelun.com

► URINE / STOOL CONTAINER SERIES



Urine Cup/Container

- Item code: CMUC
- Specification: 30ml, 40ml, 60ml, 100ml, 120ml
- Sterile or Non-sterile
- With or without label



Stool Cup/Container with Spoon

- Item code: CMSC-1
- Specification: 30ml, 40ml, 60ml, 100ml, 120ml
- Sterile or Non-sterile
- With or without label




Specimen Container

- Item code: CMSC-2
- Specification: 40ml, 60ml, 90ml
- Sterile or Non-sterile
- With or without label



Vacuum Urine Collector - Urine Collection Cup

- Item code: CMUB-3
- Size: 100ml, 120ml
- Secure collection of urine for ECBU with the vacuum tube



LeLun, Life Guard!

MASK SERIES



Nebulizer Mask

Item code: CMNM

Features:

- Made of medical PVC
- Without nose clip, face-fit shape
- Medium concentration
- Flexible connector, 360° rotation, more convenient for use
- Available with anti-crush 2m tubing
- 6cc or 10cc nebulizer cup



Oxygen Mask

Item code: CMOM

Features:

- Human body imitation design, makes a good fit with the chin, and a good air-tightness
- Transparent and clear medical grade PVC for patient comfort and visual assessment
- Adjustble nose clip which can be optional assures comfortable fit
- Available with oxygen mask without nose clip

. 60

- Latex Surgical Glove

 Item code: CMGS
- Material: Nature latex
- Specification: 6, 6.5, 7, 7.5, 8, 8.5
- Powdered and powder free
- Package: One pair/pouch, 50 pairs per box

Latex Examination Glove

- Item code: CMGE
- Material: Nature latex
- Specification: S, M, L, XL
- Grade: AQL1.5, AQL2.5, AQL4.0
- Powdered and powder freePackage: 100 pcs per box

Nitrile Examination Glove

- Item code: CMNGE
- Material: Nitrile
- Specification: S, M, L, XL



Non-rebreather Oxygen Mask (with Reservoir Bag)

Item code: CMRM

Features:

- Mask with elastic strap
- Adjustable nose clip
- With 2m tubing
- Size: XS, S, M, L, L3, XL
- Bag: 1000ml or 600ml

Nasal Oxygen Cannula

Item code: CMNOC Features:

- Size: Adult, Child, Infant, NeonateFree movement with the tube, can
- work with 15kpa. of airflow

 Available with standard or universal connector, and different lengths of oxygen tubing
- Double hole nare



Oxygen Connection Tube

Item code: CMOCT

Features:

- Free movement with the tube, can work with 15kpa. of airflow
- Match with 5mm star lumen tube, anti-crush
- Being made of all new material, enables it to be soft, smooth and elastic, without twist or knot
- Available with standard or universal connector, and different lengths



HERNIA MESH SERIES

MEDICAL GLOVE SERIES





Item code: CMHM

Features:

- Using medical non-toxic polypropylene mesh
- The surgical trauma is small, only local anesthesia is needed, and the operation is convenient
- After the mesh plug is placed in the defect area of the groin, it can effectively disperse the abdominal pressure
- It can be cut appropriately according to the size of the hernia defect
- It can quickly promote the response of fibroblasts after surgery, and grow into the gap of the mesh to form a solid fibrous tissue wall
- Can be fixed in the human body for a long time

15 / www.lelun.com



COMPONENTS

Disposable Safety Hypodermic Needle Size:

16G/17G/18G/19G/20G/21G/22G/23G/24G/ 25G/26G/27G

- Features:
- 1. Safety sheath: Safety mechanism ensures the
- 2. Standard hub: Fit with standard luer connectors, ensure qualified medicine residues
- 3. Four locking points: Totally locked, ensure safety





Drip Chamber

- Length: 48mm / 54mm
- Size: 20 drips / 1ml; 60 drips / 1ml



Combi Stopper

- Fit for both male and female luer
- Individual blister package



Syringe Filter

- Material: Hydrophilic PVDF
- Size: 0.22μm / 0.45μm
- Oiameter: 13mm / 25mm
- Ocnnector: Luer tip



Latex Connector

Used with the infusion set





Y Connector

Material: PVC



Flow Regulator with Horn Material: PE



Adaptor

Material: PE

CUSTOMERS FROM ALL OVER THE WORLD

















