
Ear Microsurgery


Malleus Nippers, Instrument Guard





CATALOG KARL STORZ








222800 HOUSE-DIETER Malleus Nipper, dismantling, upbiting, working length 8 cm

222900 Same, downbiting

222801 FISCH-DIETER Malleus Nipper, dismantling, e

FISCH-DIETER **Malleus Nipper, dismantling,** extra delicate, upbiting, width 0.8 mm, working length 8 cm



223000 **Instrument Guard,** Teflon, **autoclavable,** for micro forceps and micro scissors with working length 8 cm, length 11.5 cm


Ear Microsurgery

Scissors



Scissors, working lengths 8 cm and 9 cm








222401	WULLSTEIN Scissors , delicate, straight, with cleaning connector, working length 9 cm
222402	Same, curved upwards
222403	Same, curved to left
222404	Same, curved to right
222405	WULLSTEIN Scissors , extra delicate, straight, with cleaning connector, working length 8 cm

Scissors and Micro Scissors, working length 8 cm, dismantling





	222500	BELLUCCI Scissors, dismantling, cutting length 8 mm, working length 8 cm
	222501	FISCH-BELLUCCI Scissors, dismantling, extra fine blade, cutting length 7 mm, working length 8 cm
	222600	BELLUCCI Scissors, dismantling, delicate standard model, cutting length 8 mm, working length 8 cm
7	222601	FISCH-BELLUCCI Micro Ear Scissors, dismantling, beaked scissors blades, cutting length 7 mm, working length 8 cm
	222602	HOUSE-BELLUCCI Scissors, dismantling, extra delicate, working length 8 cm
	222604 R	BELLUCCI Scissors, dismantling, delicate, curved to right, rounded, working length 8 cm
	222605 L	Same, curved to left


(-)

2 55

Sleeve Gastrectomy 59

BERCI Fascial Closure Instrument


for subcutaneous fascial closure



Palpation Probe

Size 5 mm

Operating instruments, **length 43**, for use with laparoscopes with inbuilt working channel and **trocars size 6 mm**



Accessories

Unipolar and Bipolar High Frequency Cords



Unipolar High Frequency Cords



Unipolar High Frequency Cord, with 4 mm plug, length 300 cm, for models KARL STORZ, Erbe type T, older models and Ellman



Bipolar High Frequency Cords

26176LV



Bipolar High Frequency Cord, length 300 cm, for AUTOCON® II 400 SCB system (112, 114, 116, 122, 125), AUTOCON® II 200, AUTOCON® II 80 and Valleylab coagulators

Standard Instruments for Enucleation

HOPKINS® Telescopes






27005BA HOPKINS® Forward-Oblique Telescope 30°,

enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated,

color code: red



27005FA

HOPKINS® Telescope 12°, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated,

color code: black

Resectoscope Sheaths



27050SCK

Resectoscope Sheath, quick-release lock, including connecting tube for in- and outflow, 26 Fr., oblique beak, rotating inner sheath with ceramic insulation, color code: yellow

27050SLK

Resectoscope Sheath, including connecting tube for inand outflow, 26 Fr., oblique beak, rotating inner sheath with ceramic insulation, color code: yellow

Obturator



27040OC

Standard Obturator, for 24/26 Fr. sheaths

27050BK


SCHMIEDT **Visual Obturator**, for 24/26 Fr. sheaths, 27040BK/BO/SD/SL, 27240BO, 27241BK/BO, 27042B, 27242BZ, 27050SL

Products & Solutions Customer Service Online Catalog About us

Home page > Online Catalog - Human Medicine

27050SC





Product group / field of application

Product details

Resectoscope sheath with rotating inner sheath and click mechanism for continuous irrigation and suction

- · Sheath can be connected in any position
- · Rotating inner sheath
- · For better service: Exchangeable inner sheath

Type of usage

Technical details

	Les
System size	26 Fr.
Color code	Yellow
Flow property	Continuous flow
LUER-Lock	Rotating
Sheath end	Oblique
Connection mechanism	Click mechanism
Sheath rotation	Rotatable inner sheath
Sterile	No

Reusable

Documents & Films

Resectoscope Sheaths with Rotating Inner Sheaths and Click Mechanism


for continuous irrigation and suction



Special Features:

- Sheath can be connected in any position
- Rotating inner sheath
- For better service: Exchangeable inner sheath

 Sheaths can be combined with Working Elements 27040 and 27050



27050SCK

Resectoscope Sheath, including tube for in- and outflow, 26 Fr., oblique beak, rotating inner sheath with ceramic insulation, quick-release lock color code: yellow

27050CA

Inner Sheath, rotating, with ceramic insulation, for use with resectoscope outer sheath of 27050SC

27050SDK

Resectoscope Sheath, including tube for in- and outflow, 28 Fr., oblique beak, rotating inner sheath with ceramic insulation, quick-release lock color code: black

27050CB

Inner Sheath, rotating, with ceramic insulation, for use with resectoscope outer sheath of 27050SD

The resectoscope sheaths can be used with unipolar and bipolar working elements.


le resectoscope sheaths can be used with unipolar and bipolar working elemen

6

HOPKINS® Telescopes







7-22