

## OUR COMMITMENT:

**Hangzhou Medcare Medical Products Co.,Ltd** ,one-step solution for medical products Purchase, With an experienced team for medical field for more than 5 years,focused on supply for Disposable plastic blood bags. Our manufactory is the senior enterprise in China which was established in year 1986,and the blood bags have been sold to blood banks and donation Centres all over the domestic market such as Shanghai blood banks center, also been exported to European, African and Asia markets etc.

All manufacturing proccess are quality controlled and maintained to the highest standars.Blood management systems include all solutions such as CPDA-1,CPD+SAGM allowing platelet storage time of 5-7days or filtration system etc.

Quality management system including CE and ISO9001,ISO13485,and the anticoagulant Contained in the blood bag has been GMP certificated by SFDA, also have Apostilled Free Sale Certificate etc to meet different coutries registration.

In addition, we also export other consumables, including blood transfusion set, blood vacuum tubes, syringes, plaster of bandages,surgical and examination gloves and Dengue kits and x-ray film etc to meet the customer's demands.

We believe good quality can be long time business,welcome to your inquiries!



## HANGZHOU MEDCARE MEDICAL PRODUCTS CO.LTD

📍 Add:RM 203 Buliding no#4 East Software Park,  
#90Wensan Road, Hangzhou ,China

☎ Tel: +86 13957147193

✉ Email: Lisa@hzmedcare.com.cn

🌐 www.hzmedcare.com.cn



## Disposable plastic Blood bag

◆ 1.Single Blood bag

Single bag is used for collection of whole Blood with anticoagulant CPDA-1 solutions USP.

Availabled two styles:Tubular film(Style A)& Flat film (Style B)



Style A

Style B

Size and Specification:

| Product Code      | Anticoagulant Additive  |          | Bag capacity | Needle Gauge Size | Packaging    |
|-------------------|---|----------|--------------|-------------------|--------------|
|                   | Type  | Quantity |              |                   | Pouch/Carton |
| MS250             | CPDA-1  | 35ml     | 250ml        | 16G/17G           | 5/100        |
| MS350             |   | 49ml     | 350ml        |                   | 5/100        |
| MS450             |   | 63ml     | 450ml        |                   | 5/100        |
| MS500             |   | 70ml     | 500ml        |                   | 5/100        |
| OptionalFeatures: | Pre-donation sampling device ,needle protector,customized accessories |          |              |                   |              |

◆ 2.Double Blood bag

Double bag is used for separation of two components from whole blood, it includes one primary bag with anticoagulant CPDA-1 Solutions USP and one empty satellited bag.

Availabled two styles:Tubular film(Style A) & Flat film (Style B)



Style A

Style B

Size and Specification:

| Product Code      | Anticoagulant Additive  |          | Bag capacity | Needle Gauge Size | Packaging    |
|-------------------|---|----------|--------------|-------------------|--------------|
|                   | Type  | Quantity |              |                   | Pouch/Carton |
| MD250             | CPDA-1  | 35ml     | 250+200ml    | 16G/17G           | 4/100        |
| MD350             |   | 49ml     | 350+300ml    |                   | 4/100        |
| MD450             |   | 63ml     | 450+400ml    |                   | 4/80         |
| MD500             |   | 70ml     | 500+400ml    |                   | 4/80         |
| OptionalFeatures: | Pre-donation sampling device ,needle protector,customized accessories |          |              |                   |              |

◆ 3.Triple Blood Bag

Triple blood bag is used for the separation of three components from whole bag, including one primary bag With Anticoagulant CPDA-1 solution BP and two empty satellite bags.

The triple SAGM blood bag, includes one primary bag with anticoagulant CPD solution UPS, SAGM solution in another bag and one empty satellite bag. only suggest flat film style B.

Size and Specification:

| Product Code      | Anticoagulant Additive  |          | Bag capacity | Needle Gauge Size | Packaging    |
|-------------------|---|----------|--------------|-------------------|--------------|
|                   | Type  | Quantity |              |                   | Pouch/Carton |
| MT250             | CPDA-1  | 35ml     | 250+200*2ml  | 16G/17G           | 4/80         |
| MT350             |   | 49ml     | 350+300*2ml  |                   | 4/80         |
| MT450             |   | 63ml     | 450+400*2ml  |                   | 3/72         |
| MT500             |   | 70ml     | 500+400*2ml  |                   | 3/72         |
| MT250S            | CPD-SAGM  | 35/50ml  | 250+200*2ml  | 16G/17G           | 3/60         |
| MT350S            |   | 49/75ml  | 350+300*2ml  |                   | 3/60         |
| MT450S            |   | 63/100ml | 450+400*2ml  |                   | 3/54         |
| MT500S            |   | 70/100ml | 500+400*2ml  |                   | 3/54         |
| Optional Features | Pre-donation sampling device,needle protector, 5-day platelet bag,filtration system etc custmized accessories |          |              |                   |              |



Style A

Style B

◆ 4.Quadruple Blood bag /Paediatric Blood bag


Quadruple blood bag is used for the separation of four components from whole human blood , The quadruple system includes one primary bag with anticoagulant CPDA-1 solution USP and three empty satellite bags.

The Quadruple SAGM blood bag includes one primary bag with anticoagulant CPD solution USP,SAGM Solution in another bag and two.empty satellite bags

Paediatric Blood bag: 350ml or 450ml one primary bag with three 100ml or 150ml satellite bags

Size and Specification:

| Product Code      | Anticoagulant Additive   |          | Bag capacity | Needle Gauge Size | Packaging    |
|-------------------|--|----------|--------------|-------------------|--------------|
|                   | Type   | Quantity |              |                   | Pouch/Carton |
| MQ250             | CPDA-1   | 35ml     | 250+200*3ml  | 16G/17G           | 3/66         |
| MQ350             |  | 49ml     | 350+300*3ml  |                   | 3/66         |
| MQ450             |  | 63ml     | 450+400*3ml  |                   | 3/60         |
| MQ500             |  | 70ml     | 500+400*3ml  |                   | 3/60         |
| MQ250S            | CPD-SAGM   | 35/50ml  | 250+200*3ml  | 16G/17G           | 3/60         |
| MQ350S            |  | 49/75ml  | 350+300*3ml  |                   | 3/60         |
| MQ450S            |  | 63/100ml | 450+400*3ml  |                   | 3/54         |
| MQ500S            |  | 70/100ml | 500+400*3ml  |                   | 3/54         |
| Optional Features | Pre-donation sampling device,needle protector, 5-day platelet bag,filtration system etc custmized accessories. |          |              |                   |              |



Style B

MHF01

◆ 5.Disposable Blood Transfusion Set

95mm drip with filter inside ;Roller clamp; Latex flash ball with needle connect to luer lock 150cm tubing in total length  
EO sterilization