

S130D PM KIT-1 (Every 600 cycles)

No.	Part	Picture	Q'ty	Spec.	Cycle	Year	
	PM06-0010-00 (Option 1)						
	PM06-0020-00 (Option 2)						
1-1	Rotary Pump Oil (Laybold)		1	1 Liter	600	0.5	
1-2	Rotary Pump Oil (WooSung)		3	1 Liter	600	0.5	
2	Activated Carbon Filter		1		600	0.5	
3-1	O-ring (Inlet&Outlet) (Laybold)		1		600	0.5	
3-2	O-ring (Inlet&Outlet) (WooSung)		1		600	0.5	
4	Oil Paper & Glove		1		600	0.5	
5	Oil container		1	3 Liter	-	-	
6	Funnel		1		-	-	
7	Garbage Vinyl		1		-	-	

※ Remark

Replace PM KIT-1 every 6 month or 600 cycles, whichever comes first.

※ Cleaning

Injection port hole cleaning should be performed every 600 cycles.

S130D PM KIT-2 (Every 1,200 cycles)

Pri No.	Part	Picture	Q'ty	Spec.	Cycle	Year	
	PM06-0030-00 (Option 1)						
	PM06-0040-00 (Option 2)						
	PM06-0060-00 (Option 3)						
1-1	Option 1 Oil Mist Trap (Laybold)		1		1,200	1.0	
1-2	Option 2 Oil Mist Trap (WooSung)		1	SUS	1,200	1.0	
1-3	Option 3 Oil Mist Trap-Inside Filter Set (WooSung)		1	Set	1,200	1.0	
2	Heat-registing O-ring		1		1,200	1.0	
3	NW25 Clamp		1		-	-	

※ Remark

Replace PM KIT-2 every 1 year or 1,200 cycles, whichever comes first.

※ Cleaning

Injection ball inside of SFS Nozzle cleaning should be performed every 1,200 cycles.

S130D PM KIT-3 (Every 2,400 cycles)

No.	Part	Picture	Q'ty	Spec.	Cycle	Year	
	PM06-0050-00						
1	Catalyst Filter		1		2,400	2.0	
2	Injection Valve		2	Shintec	2,400	2.0	
3	HEPA Filter		1		2,400	2.0	
							

※ Remark

Replace PM KIT-3 every 2 years or 2,400 cycles, whichever comes first.