


Table Top Low Speed Centrifuge



Introduction:

BKC-TL5B is used for separation of serum and plasma, precipitation of protein or urine sediment examination.

Features:

- * It can store 100 programs, automatic memory form the lasting running programs.
- * Adjustable acceleration levels: 20; deceleration levels: 20.
- * All-metal structure multi-layer explosion-proof design, strong and durable.
- * Embedded microprocessor control, DC brushless motor drive, maintenance free, smooth operation.
- * 5 inch touch LCD screen, real-time display all running parameters, the operation interface is intuitive and simple.
- * Support small language customization,
- * Mechanical and electrical integration door lock, easy to use, safe and reliable.
- * Silicone integral sealing ring to avoid aerosol overflow.
- * With self-test, block, overspeed, imbalance, over temperature detection and alarm functions, excellent shock absorption device.
- * Store multiple groups of user-defined programs, boot for the last program.

Technical Parameters:

Model	BKC-TL5B
Max.Speed	5000rpm
Max.RCF	3894×g
Max.Capacity	4*250ml
Speed Precision	±10rpm
Time Range	1s~99h59min/Continuous/Short-span
Power Supply	AC110/220V 50/60Hz
Noise	≤65dB(A)
External Size(L*W*H)	430*565*410mm
Package Size(L*W*H)	530*665*570mm
Net/Gross Weight	43/52kg

Rotor:

Type	NO.	Capacity	Max.Speed(rpm)	Max.RCF(×g)	Adapter
Swing Rotor	5M00024	48*2/3/4/5/7ml	4000	2952	/
	5M00025	48*5/7ml		2952	2/3/4ml vacuum blood collection tube
	5M00026	72*2/3/4/5/7ml		2808	/
	5M00029	16*15/10ml		2957	2/3/4ml vacuum blood collection tube,
	5M00030	24*15/10ml		2952	5/7ml vacuum blood collection tube
	5M00031	32*15/10ml		2952	5ml, 7ml,15ml(▲)
	5M00032	4*50ml	4500	3747	50ml(▲), 20ml, 15ml, 10ml, 5ml, 1.5ml
	5M00033	8*50ml	4000	2952	50ml(▲), 20ml, 15ml, 10ml, 5ml, 1.5ml
	5M00034	4*100ml	4500	3894	2/3/4ml vacuum blood collection tube,
	5M00036	8*100ml	4000	3077	5/7ml vacuum blood collection tube
	5M00038	4*250ml	4000	2780	48*1.5/2ml, 40*2/3/4/5/7ml vacuum blood collection tube, 24*10ml(▲),28*10ml, 20*15ml,4*50ml(▲), 8*50ml(●), 4*100ml
	5M00058	4*300ml	4000	3220	/
Microplate Rotor	5M00045	2*2*96well	4000	2361	/
Angle Rotor	5M00072	12*20ml	5000	3516	15ml, 10ml, 5ml, 1.5ml
	5M00073	12*15ml	5000	3354	1.5ml
	5M00074	8*15ml	5000	2655	1.5ml
	5M00075	24*10ml	5000	3354	1.5ml, 2/3/4ml vacuum blood collection tube
	5M00076	18*10ml	5000	3354	
	5M00077	12*10ml	5000	3354	
	5M00084	6*50ml	5000	2711	50ml(▲), 20ml, 15ml, 10ml, 5ml, 1.5ml