



# Adhesive fabric tape SENSIplast

TD-06-I.2.c-5.5b

edition I

date: 2021-05-26

amendment of: 2021-12-22

**MANUFACTURER:** ZARYS International Group sp. z o.o. sp.k.**ADDRESS:** ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com**CLASSIFICATION**

wyrób medyczny klasy I zgodny z Rozporządzeniem PE i Rady (UE) 2017/745 (MDR)

**INTENDED PURPOSE**

- device intended for stabilisation, protection and fixation of drains, catheters, medical tubes not requiring specialised fixation at their place of application as well as primary dressings (contact), non-adhesive, applied directly to wounds and skin damages

**FEATURES**

- made of cotton fabric
- white colour
- wound on a spool
- covered with an acrylic adhesive applied by the continuant method
- can be divided lengthwise and crosswise without the use of scissors
- serrated edge finishing
- good air and water vapour permeability
- good adhesion
- thin, flexible, conforms to the body contours, guarantees freedom of movement
- hypoallergenic
- single-use
- non-sterile
- shelf life: 3 years



CATALOGUE NUMBER	DIMENSIONS	INTERMEDIATE PACKAGING (SALES) BOX	BULK PACKAGING CARTON
STA12505W	1,25 cm x 5 m	24 x 1 pc.	30 x 24 pcs
STA12509W	1,25 cm x 9,14 m	24 x 1 pc.	30 x 24 pcs
STA25005W	2,5 cm x 5 m	12 x 1 pc.	30 x 12 pcs
STA25009W	2,5 cm x 9,14 m	12 x 1 pc.	30 x 12 pcs
STA50005W	5 cm x 5 m	6 x 1 pc.	30 x 6 pcs
STA50009W	5 cm x 9,14 m	6 x 1 pc.	30 x 6 pcs



elastopor STERIL  
non-woven dressing with absorbent pad,  
self-adhesive, sterile

TD-09

att.no.3 pt.2 subpt.2.3

edition I

date: 2021-10-29

change of day: 2021-12-22

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

**INDICATION:**

- treating wounds with a small and medium amount of exudate
- especially recommended for covering surgical wounds characterized by straight cuts and the application of fewer or more skin sutures in a straight line

**FEATURES:**

- self-adhesive dressing with absorbent pad
- made of hydrophobic non-woven material
- micropores in the non-woven material structure provide an adequate rate of moisture vapour permeability so correct gas exchange between the dressing and skin is guaranteed
- centrally placed absorbent pad - absorbs medium wound exudate
- absorbent pad coated with polyethylene mesh protects against external factors' influence and does not stick to the wound as well as facilitates removal of the dressing
- hypoallergenic
- elastic - comfortable for the patient, adjusts to the body shape and guarantees freedom of movement
- easy to attach to the skin
- rounded edges prevent accidental removal of the dressing
- covered with hypoallergenic acrylic glue without the content of zinc oxide, rubber, latex, evenly applied over the entire surface, with very good moisture resistance ensuring good stabilization of the dressing at the application site
- silicone paper backing facilitates precise, painless and efficient application
- EO sterilized



**SIZE AND REFERENCE**

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-PAPER	INTERMEDIATE PACKAGING BOX	BULK PACKAGING CARTON
801001E	5 cm x 7,2 cm	1 pc.	100 pcs	24 x 100 pcs
801012E	6 cm x 10 cm	1 pc.	50 pcs	45 x 50 pcs
801002E	8 cm x 10 cm	1 pc.	30 pcs	36 x 30 pcs
801013E	10 cm x 10 cm	1 pc.	30 pcs	24 x 30 pcs
801021E	10 cm x 12 cm	1 pc.	30 pcs	30 x 30 pcs
801014E	10 cm x 15 cm	1 pc.	30 pcs	24 x 30 pcs
801003E	6 cm x 15 cm	1 pc.	50 pcs	24 x 50 pcs
801004E	8 cm x 15 cm	1 pc.	30 pcs	24 x 30 pcs
801005E	8 cm x 20 cm	1 pc.	30 pcs	18 x 30 pcs
801015E	10 cm x 20 cm	1 pc.	25 pcs	18 x 25 pcs
801016E	10 cm x 25 cm	1 pc.	25 pcs	18 x 25 pcs
801017E	10 cm x 30 cm	1 pc.	25 pcs	12 x 25 pcs
801018E	10 cm x 35 cm	1 pc.	25 pcs	12 x 25 pcs