



GEISTER®

ValveGate™ PRO



Efficient surgery. Less invasive.
Optimum results.





ValveGate™ TC Graspers 2.7 mm shaft

- A 38-7935 10° angled
- B 38-7936 25° angled
- 38-7934 straight

ValveGate™ PRO Graspers

- C 38-7824 Resano, strong, jaw for valves, straight
- D 38-7929 Russian, jaw 6.8x5.3 mm
- E 38-7567 Lymph Node Grasper angled up, jaw 7 mm ovale, atraumatic

ValveGate™ PRO Gentle Graspers

- F 38-7816 extra wide jaw, straight, TC
- G 38-7826 soft extra wide jaw, straight, 2x3 teeth, 4.8x14 mm, with extra fine teeth

ValveGate™ PRO Graspers with DeBakey Jaws, 1x2 teeth

- H 38-7820 narrow, straight, 1.5x11 mm
- I 38-7821 narrow, angled upwards, 1.5x14 mm
- J 38-7822 narrow, angled sideways, 1.5x18 mm, vertically opening
- K 38-7825 wide jaw, 2.8x11 mm, straight
- L 38-7827 delicate, 1.5x11 mm, 3 mm malleable shaft
- M 38-7881 narrow, 1.5x15 mm, straight
- N 38-7917 narrow, 1.5x15 mm, slightly curved
- O 38-7918 narrow, 1.5x15 mm, angled sideways 10°
- P 38-7919 narrow, 1.5x15 mm, angled sideways 25°
- Q 38-7922 narrow, angled sideways, 1.5x18 mm, horizontally opening
- R 38-7928 long, 2.8x25 mm teeth, 40 mm jaw

5 mm shaft diameter, unless otherwise specified. Standard working length 25 cm/10" for minimally invasive valve surgery. For different length add "B", "C", "S", "L" or "XL" behind article number.

B = 12 cm/5" C = 15 cm/6" S = 17 cm/6¾" L = 30 cm/12" XL = 35 cm/14"