

TECHNICAL DATA SHEET

TYVEK REELS

PRODUCT CODE : 4ATR

RELATED STANDARDS : 4A Medical Tyvek Reels are in compliance with the provisions of the Medical Device Regulation

2017/745 for medical devices.

INTENDED USE : The 4A Medical Tyvek Reels are used to sterilize temperature-sensitive medical devices. It is preferred

at Low Temperature Sterilization regarding to its porous structure. Reels meet the Medical Industry's requirement for high quality, hygiene and safety levels. 4A Medical reels are used for the packing of the Medical Devices before sterilization and it maintains the sterility of the products after sterilization.

OPERATION CONDITIONS : 4A Medical Tyvek Reels are used for gas plazma (VH2O2) sterilization. The sterilization conditions

should be determined by the end user regarding to material to be sterilized.

SPECIFICATIONS

	PROPERTIES	UNITS	TYPICAL VALUE
	DELAMINATION	N/2.54 cm	2,2
	BASIC WEIGHT	g/m²	64
	POROSITY	sec	22
	OPACITY	%	91
	THICKNESS	mils	165
TYVEK	SPENCER PUNCTURE	J/m	6195
Ę.	MULLEN BURST	kPa	1055
	HYDROSTATIC HEAD	cm H2O	145
	ELMENDORF TEAR, XD	mN	3203
	ELMENDORF TEAR, MD	mN	2980
	TENSILE STRENGTH,XD	N/2.54 cm	174
	TENSILE STRENGTH,MD	N/2.54 cm	163
	FREE FROM LEAD	D AND HEAVY METALS AND	TOXIC MATERIALS

TEHNICALDATA SHEET



TYVEK REELS

NATED FILM	PROPERTIES	UNIT	TYPICAL VALUE		
			LAM.FILM	52 ± 5	
	THICKNESS	Micron	PE	40 ± 5	
			PET	12 ± 2	
	BREAKING FACTOR	N/15mm	35		
	TEAR STRENGTH MD/CD	mN	300		
MINA	BOND STRENGHT	N/m	265		
LA	ELONGATION	%	70		
	TEMPERATURE DURABILITY	min/degrees	30 min/140C		
	FREE FROM LEAD AND HEAVY METALS AND TOXIC MATERIALS				

PRODUCT SPECIFICATION

	PROPERTIES	UNIT	TYPICAL VALUE
	SEAL STRENGTH	N/15mm	2,5
	BUBBLE TEST		None
СН	PINHOLE		None
P01	DIMENSION CONTROL	cm or mm	Desired dimensions
-	LEAKAGE TEST		None
	PEEL DIRECTION		None
	INDICATOR CONTROL		Returns the specified color

CATO	STERILIZATION METHOD VH2O2		BEFORE STERILIZATION	AFTER STERILIZATION
INDICATO R			Red	Yellow
	PACKAGING	:	4A Medical reels are packaged in carton	boxes as below:
		•	Product dimensions	Pieces in carton
			75mm X 70m	8
			100mm X 70m	8
			150mm X 70m	4
			200mm X 70m	4
			250mm X 70m	2
			300mm X 70m	2
			350mm X 70m	2
			400mm X 70m	2
			450mm X 70m	2
			500mm X 70m	2
MANI	FACTURER	:	DORT-A TIP MALZEMELERI SAN. ITH. IHR. TIC. LTD. STI.	
			Balıkhisar Mahallesi, Köyiçi Serpmeleri, No:795 Akyurt, Ankara, Turkey Tel: +90 312 363 50 52 / 53 Fax: +90 312 363 52 10	
			Web: http://www.4amedical.com E-mail: export@4amedical.com	

TECHNICAL DATA SHE

E