

#### **Summary Technical Documentation**

TD1029\_8 Revision:8 Page 1 of 1

# **DECLARATION OF CONFORMITY**

optimed Medizinische Instrumente GmbH Ferdinand-Porsche-Straße 11 76275 Ettlingen

hereby declares in its own responsibility that

### INTERVENTIONAL NEEDLES

Order codes:
Please refer to the annexed article range
of TD1029\_1.1.7, Rev. 5

is a Class IIa medical device reviewed in accordance with Annex II excluding section (IV)

0124

DEKRA Certification GmbH Handwerkstrasse 15 D-70565 Stuttgart

conforming to the relevant provisions of the Directive 93/42/EEC.

This Declaration of Conformity is valid until 26.05.2024.

The validity corresponds to the validity of the CE-certificate number 50066-16-08.

Ettlingen, M.OS. 2021

Frank Wohkittel

Head of Regulatory Affairs

noptimed (

Medizinische Instrumente GmbH Ferdinand-Porsche-Strasse 11 76275 Ettlingen · Germany Phone +49 7243/7633-0

Fax +49 7243/7633-624







#### **Summary Technical Documentation**

TD1029\_1.1.7 Revision:5 Approved on: 04.12.2020 21:59 Page 1 of 2

Person in charge: Christian Eing Last amended by: Ruth Siewe Reviewed by: Unterliegt keiner Prüfung Approved by: Zlatka Schüßler-Tankova

## 1.1.7 Complete List of the Various Configurations and Variants

The range of interventional needles offered by optimed consists of the following categories:

- Angiography Needles
  - o US-cut
- Initial Puncture Needles (with or without marking)
  - o US-cut
- Chiba Needles (with or without marking)
  - o US-cut
  - o bevel-edged

Extension	Reference	Specification	
_AGN	1201-1200	Angiography Needle; Needle, 1-part, US-cut; D: 18 G (1.2 mm), L: 70 mm	
_IPN	1403-0151	Initial Puncture Needle; Needle, 1-part, US-cut; D. 17.5 G (1.3 mm), L: 150 mm	
_IPN	1403-0200	Initial Puncture Needle; Needle, 2-part, hollow stylet: US-cut, outer cannula: blunt tapered, marking; D: 17.5 G (1.3 mm), L: 200 mm	
_CHN	1411-1150	Chiba Needle; Needle, 2-part, full stylet: bevel edged, outer cannula: bevel edged, marking; D: 22 G (0.70 mm), L: 150 mm	
_CHN	1411-1220	Chiba Needle; Needle, 2-part, full stylet: bevel edged, outer cannula: bevel edged, marking; D: 22 G (0.70 mm), L: 220 mm	
_CHN	1413-1090	Chiba Needle; Needle, 2-part, full stylet: round blunt, outer cannula: US-cut, marking; D: 22 G (0.70 mm), L: 90 mm	
_CHN	1413-1150	Chiba Needle; Needle, 2-part, full stylet: round blunt, outer cannula: US-cut, marking; D: 22 G (0.70 mm), L: 150 mm	
_CHN	1413-1220	Chiba Needle; Needle, 2-part, full stylet: round blunt, outer cannula: US-cut, marking; D: 22 G (0.70 mm), L: 220 mm	
_CHN	1413-3220	Chiba Needle; Needle, 2-part full stylet: round blunt, outer cannula: US-cut; D: 19.5 G (0.95 mm), L: 220 mm	

#### Key:

AGN	Angiography Needle	
IPN	Initial Puncture Needle	
_CHN	Chiba Needle	

The following raw material is not being sold as single item but as a part of another product (Drainage Set).

Reference	Specification
-----------	---------------







## **Summary Technical Documentation**

TD1029\_1.1.7 Revision:5 Approved on: 04.12.2020 21:59 Page 2 of 2

012209090 Chiba Needle; Needle, 2-part, full stylet: round blunt, outer cannula: US-cut; D: 19.5 G (0.95 mm), L: 90 mm



