

## BioPlast Special

**IVD** In vitro diagnostic – medical device

CND code: W01030799

### General features

Mixture of pure paraffins and other waxes, specific for the processing of all kind of samples.

Code	Size	Packaging
08-7915	2 kg	6 bags of 2 kg each

The product validity correctly stored is 5 years.

Manufacturer: **Bio-Optica Milano S.p.A.**

### Technical data

<b>Melting point</b>	<b>52 – 54°C</b>
Oil percentage (ASTM D721)	About 0,5%
Shape	Drops
Colour	White



Use destination For processing

The drop's shape allows an easy handling and promotes melting times.

BioPlast Special wax is manufactured without any toxic contaminant, also because it has a triple filtration during the productive passage that makes it pure.

The product's high quality is guaranteed by an automated productive process really effective.

### Warnings and use's method

- At temperatures upper than 63°C, the waxes lose their plastic properties.  
Check always the waxes' temperatures in processors and embedding stations.
- **It doesn't contain DMSO.**

**NOTE:** for further information regarding risks, the first help measures and precaution measures to take during the use and the storage, it's available the safety data sheet.