

Vaculine Maxi and Vaculine Maxi E Surgical Aspirator

Ref: M4E 500 0 0 - Vaculine Maxi

Ref: M4E 650 0 0 - Vaculine Maxi E

For operating room, critical care, emergency care and outpatient clinic use

- ▶ 2, 3 and 5 L single/double polycarbonate bottle combinations
- ▶ Suitable for use with canister systems
- ▶ 0-680/710 mmHg vacuum range (at sea level)
- ▶ Vacuum diversion knob to select bottle
- ▶ Vacuum adjust knob and 100mm gauge
- ▶ Output: 50L/ml Air filter
- ▶ One suction bag remains functional while other is being discarded *
- ▶ Ergonomic control panel design with 25 degree angle
- ▶ Fluid aspiration into the pump:
 - Floating valve,
 - Hydrophobic filter:
 - Tight seal lid with clear view of the filter, enabling prompt filter change
 - Hydrophobic filter allows free air-flow, but prevents fluid overflow, protecting the pump
- ▶ Steering handle
- ▶ Basket in the rear for spare suction catheters*
- ▶ Rinse solution receptacle and two separate suction line holders on the sides
- ▶ Polycarbonate bottle/canister compatible with steam sterilization (Applicable for 2 and 3 L bottles/canisters)
- ▶ Operation status indication lamp
- ▶ Silent, dry vacuum pump
- ▶ Design features to prevent airborne contamination in operating rooms and clean rooms. No requirement for periodic lubrication and maintenance
- ▶ Motor speed: 1400 rpm
- ▶ Sturdy body construct with durable polyester dye
- ▶ Mobile
- ▶ Four 360° swivel castors with 100 mm diameter, two with brakes
- ▶ Rated voltage: 230 V, 50/60 Hz
- ▶ Power consumption: 145W
- ▶ Dimensions (Vaculine Maxi): 485 x 650 x 985 mm
- ▶ Dimensions (Vaculine Maxi E): 485 x 655 x 965 mm
- ▶ Weight (Vaculine Maxi): 35 kg
- ▶ Weight (Vaculine Maxi E): 32 kg
- ▶ Packing dimensions: 630 x 630 x 1080 mm

Vaculine Maxi and Vaculine Maxi E Surgical Aspirator is ordered as the basic unit. The bottle kits configuration can be selected according to customer requirements and shall be ordered separately.



Vaculine Maxi



Vaculine Maxi

Vaculine Maxi E

* Available for Vaculine Maxi; not available for Vaculine Maxi E